# ABTHERA ADVANCE™

## OPEN ABDOMEN DRESSING

## INSTRUCTIONS FOR USE

Not for use with Instillation Therapy

**DE** - Nicht zur Verwendung mit Instillation Therapy **NL** - Niet te gebruiken bij instillatietherapie  
**FR** - Ne pas utiliser avec une thérapie d'instillation **IT** - Da non usare per la terapia a instillazione  
**ES** - No debe utilizarse con terapia de instilación **DA** - Må ikke anvendes sammen med instillationsterapi  
**SV** - Inte avsett för instillationsbehandling **PTBR** - Não deve ser usado com Terapia de Instilação  
**TR** - İnstilasyon Tedavisi ile kullanıma yönelik değildir **EL** - Δεν προορίζεται για χρήση με θεραπεία ενστάλαξης  
**FI** - Ei käytettäväksi huuhteluhoidossa **NO** - Skal ikke brukes med instilleringsbehandling

Only for use with Negative Pressure Therapy provided by INFOV.A.C.™ or V.A.C.ULTA™ Therapy Units

**DE** - ABTHERA ADVANCE™ Wundauflage für den Abdominalbereich. Nur zur Verwendung mit Unterdrucktherapie der Therapiesysteme INFOV.A.C.™ oder V.A.C.ULTA™  
**NL** - ABTHERA ADVANCE™-wondverband voor open buikwonden. Alleen voor gebruik met negatieve-druktherapie geleverd door INFOV.A.C.™ of V.A.C.ULTA™ behandelingsunits **FR**  
- Pansement pour abdomen ouvert ABTHERA ADVANCE™. À utiliser uniquement avec les unités de thérapie par pression négative INFOV.A.C.™ ou V.A.C.ULTA™  
**IT** - ABTHERA ADVANCE™ medicazione per addome aperto. Esclusivamente per l'uso con la terapia a pressione negativa erogata dalle unità terapeutiche INFOV.A.C.™ o V.A.C.ULTA™  
**ES** - Apósito para abdomen abierto ABTHERA ADVANCE™. Para uso exclusivo con terapia de presión negativa suministrada mediante unidades de terapia INFOV.A.C.™ o V.A.C.ULTA™ **DA** - ABTHERA ADVANCE™ Forbinding til det åbne abdomen. Kun til brug med undertryksterapi fra INFOV.A.C.™ eller V.A.C.ULTA™ Terapienheder  
**SV** - ABTHERA ADVANCE™-förband för öppen buk. Får endast användas till behandling med negativt tryck från INFOV.A.C.™- eller V.A.C.ULTA™-behandlingsenheter  
**PTBR** - Curativo para Abdome Aberto ABTHERA ADVANCE™. Somente para uso com Terapia com Pressão Negativa aplicadas pelas Unidades de Terapia INFOV.A.C.™ ou V.A.C.ULTA™  
**TR** - ABTHERA ADVANCE™ Açık Abdomen Pansumanı. Yalnızca INFOV.A.C.™ veya V.A.C.ULTA™ Terapi Üniteleri tarafından sağlanan Negatif Basıncılı Terapi ile kullanıma yöneliktir  
**EL** - Επιδεσμός ανοικτής κοιλιακής χώρας ABTHERA ADVANCE™. Μόνο για χρήση με τα συστήματα θεραπείας αρνητικής πίεσης INFOV.A.C.™ και V.A.C.ULTA™  
**FI** - Avoimen vatsaontelon ABTHERA ADVANCE™-sidos. Käytettäväksi vain INFOV.A.C.™- tai V.A.C.ULTA™-hoitoyksiköllä annettavan alipaineimuhoidon yhteydessä  
**NO** - ABTHERA ADVANCE™-forbindingen for åpent abdomen. Bare for bruk i behandling med undertrykk med behandlingsapparatene INFOV.A.C.™ eller V.A.C.ULTA™

# ABTHERA ADVANCE™ OPEN ABDOMEN DRESSING

## INSTRUCTIONS FOR USE

Only for use with Negative Pressure Therapy provided by INFOV.A.C.™ or V.A.C.ULTA™ Therapy Units

### PRODUCT DESCRIPTION

The ABTHERA ADVANCE™ Open Abdomen Dressing, when used with Negative Pressure Therapy provided by the INFOV.A.C.™ and V.A.C.ULTA™ Therapy Units provides an active temporary abdominal closure system, designed to remove fluids from the abdominal cavity and draw wound edges together, helping to achieve primary fascial closure while protecting abdominal contents from external contaminants.

### SAFETY INFORMATION

**IMPORTANT:** As with any prescription medical device, failure to consult a physician and carefully read and follow all therapy unit and dressing instructions and safety information prior to each use may lead to improper product performance and the potential for serious or fatal injury. Do not adjust therapy unit settings or perform therapy application without directions from / or supervision by the clinical caregiver.

All disposable components of the ABTHERA ADVANCE™ Open Abdomen Dressing are for single use only. Re-use of disposable components may result in wound contamination, infection and / or failure of the wound to heal.

### INDICATIONS FOR USE

The ABTHERA ADVANCE™ Open Abdomen Dressing is indicated for temporary bridging of abdominal wall openings where primary closure is not possible and / or repeat abdominal entries are necessary. The intended use of this dressing is in open abdominal wounds with exposed viscera including, but not limited to, abdominal compartment syndrome. The intended care setting is a closely monitored area within the acute care hospital, such as the ICU. The abdominal dressing will most often be applied in the operating theater.

### CONTRAINDICATIONS

- **Never** place exposed foam material directly in contact with exposed bowel, organs, blood vessels or nerves. Protect vital structures with the Visceral Protective Layer **at all times** during therapy.
- Patients with open abdominal wounds containing non-enteric unexplored fistulas should not be treated with the ABTHERA ADVANCE™ Open Abdomen Dressing.

Management of the open abdomen has been documented in case reports and consensus panel literature. Please refer to the **References List** section of this document.

### WARNINGS

**Not for use with Instillation Therapy:** Although it is accepted medical practice to flush a contaminated open abdominal cavity with saline or other medical solutions, the ABTHERA ADVANCE™ Open Abdomen Dressing was not designed for this purpose, and KCI has no studies to support its safe and effective use with instillation therapy. Potential risks of instillation into the open abdomen include:

- Instillation of fluid in the abdomen without sufficient fluid recovery may lead to abdominal compartment syndrome.
- Instillation of fluids in the abdomen that are untested for safety and efficacy with this application could lead to severe hollow viscus and solid organ damage.
- Instillation of unwarmed fluid in large quantities may lead to hypothermia.

**Only Use the SENSAT.R.A.C.™ Pad:** Substitution with any other tubing, alteration of the SENSAT.R.A.C.™ Pad or breach of the prescribed SENSAT.R.A.C.™ Pad application for the purpose of instilling fluids into the open abdomen is not recommended under any circumstance. This may lead to loss of system efficacy or harm to the patient.

**Bleeding:** Patients with abdominal wounds must be closely monitored for bleeding as these wounds may contain hidden blood vessels which may not be readily apparent. If sudden or increased bleeding is observed in the dressing, tubing or canister, immediately discontinue Negative Pressure Therapy, take appropriate measures to stop bleeding, and contact the physician. Negative Pressure Therapy is not designed to prevent, minimize or stop bleeding.

**Hemostasis must be achieved prior to dressing placement.**

The following conditions may increase the risk of potentially fatal bleeding.

- Suturing and / or anastomoses
- Trauma
- Radiation

- Inadequate wound hemostasis
- Non-sutured hemostatic agents (for example, bone wax, absorbable gelatin sponge or spray wound sealant) applied in the abdomen may, if disrupted, increase the risk of bleeding. Protect against dislodging such agents.
- Infection in the abdominal wound may weaken visceral organs and associated vasculature, which may increase susceptibility to bleeding.
- Use of anticoagulants or platelet aggregation inhibitors.
- Bone fragments or sharp edges could puncture vessels or abdominal organs. Beware of possible shifting in the relative position of tissues, vessels or organs within the abdominal wound that might increase the possibility of contact with sharp edges.

**Intra-abdominal Pressure Monitoring:** Laparotomy with the placement of any temporary abdominal closure **does not** eliminate the possibility of elevation in intra-abdominal pressure (IAP). When using Negative Pressure Therapy, IAP monitoring (for clinical or diagnostic signs and symptoms of elevated IAP) should continue as indicated by patient condition and in accordance with institutional clinical practice or guidelines. If intra-abdominal hypertension (IAH) or abdominal compartment syndrome (ACS) is observed or suspected, note intra-abdominal pressures and turn off power to the Negative Pressure Therapy Unit, discontinuing negative pressure. After full expansion of the perforated foam, obtain a new intra-abdominal pressure measurement. If IAH / ACS persists without negative pressure, discontinue the use of Negative Pressure Therapy and address the underlying condition as medically indicated.

**Use of Visceral Protective Layer:** When using Negative Pressure Therapy, ensure that the Visceral Protective Layer completely covers all exposed viscera and completely separates the viscera from contact with the abdominal wall. Place the Visceral Protective Layer over the omentum or exposed internal organs, and carefully tuck it between the abdominal wall and internal organs, making sure the Visceral Protective Layer completely separates the abdominal wall from the internal organs.

**Adhesions and Fistula Development:** Formation of adhesions of the viscera to the abdominal wall may reduce the likelihood of fascial reapproximation and increase the risk of fistula development which is a common complication in patients with exposed viscera.

**Infection:** Infected abdominal wounds should be monitored closely and may require more frequent dressing changes than non-infected wounds, dependent upon factors such as patient condition, wound condition and treatment goals. Refer to dressing application instructions for details regarding dressing change frequency.

**Dressing Placement:** Always use a dressing from a sterile package that has not been opened or damaged. Do not force any dressing component into the wound, as this may damage underlying tissue.

**Dressing Removal:** The dressing components are not bioabsorbable. Always remove all dressing components from the abdomen at every dressing change.

**Keep Negative Pressure On:** Never leave the dressing in place without active negative pressure for more than two hours. If negative pressure is off for more than two hours, change dressing as shown in the dressing application instructions. Either apply a new dressing from an unopened sterile package and restart negative pressure, or apply an alternative dressing.

**Defibrillation:** Remove adhesive drape from area of defibrillation to prevent inhibition of electrical energy transmission.

**Acrylic Adhesive:** The drape has an acrylic adhesive coating, which may present a risk of adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives. If a patient has a known allergy or hypersensitivity to such adhesives, do not use the dressing. If any signs of allergic reaction or hypersensitivity develop, such as redness, swelling, rash, urticaria or significant pruritus, discontinue use and ensure appropriate emergency medical treatment. If bronchospasm or more serious signs of allergic reaction appear, remove dressing and ensure appropriate emergency medical intervention as indicated.

**Magnetic Resonance Imaging (MRI) – Therapy Unit:** The Negative Pressure Therapy Unit is MR unsafe. Do not take the device into the MR environment.

**Magnetic Resonance Imaging (MRI) – ABTHERA ADVANCE™ Open Abdomen Dressing:** The dressing can remain on the patient with minimal risk in an MR environment, assuming that use of Negative Pressure Therapy is not interrupted for more than two hours; please refer to **Keep Negative Pressure On** section.

**Hyperbaric Oxygen Therapy (HBO):** Do not take the negative pressure therapy unit into a hyperbaric oxygen chamber. The negative pressure therapy unit is not designed for this environment, and **should be considered a fire hazard**. After disconnecting the negative pressure therapy unit, either (i) replace the dressing with another HBO compatible material during the hyperbaric treatment, or (ii) cover the unclamped end of the SENSAT.R.A.C.™ Pad Tubing with dry gauze. For HBO therapy, the tubing must not be clamped. Never leave a dressing in place without active negative pressure for more than two hours (refer to **Keep Negative Pressure On** section).

**Application Setting:** Dressing applications and changes should be performed under strict sterile conditions in the operating theater. If dressing change is performed outside the operating theater, it must be performed in an environment equipped to address the onset of critical complications (refer to **WARNINGS** section) and where strict aseptic technique can be utilized.

### PRECAUTIONS

**Standard Precautions:** To reduce the risk of transmission of bloodborne pathogens, apply standard precautions for infection control with all patients, per institutional protocol, regardless of their diagnosis or presumed infection status. In addition to gloves, use gown and goggles if exposure to body fluids is likely.

**Intra-abdominal Packing:** When using intra-abdominal packing with Negative Pressure Therapy, packing material may be drier than anticipated. Evaluate packing material prior to removal and rehydrate if necessary to prevent adherence or damage to adjacent structures.

**Monitor Fluid Output:** The dressing is designed to efficiently remove fluid from the abdominal compartment and to evenly distribute negative pressure. When treating patients with Negative Pressure Therapy, the volume of exudate in the canister and tubing should be frequently examined.

**Patient Size and Weight:** The size and weight of the patient should be considered when prescribing Negative Pressure Therapy. Initial lower negative pressure should be considered for certain small or elderly patients who are at risk of fluid depletion or dehydration. Monitor fluid output including the volume of exudate in both the tubing and canister. This therapy has the potential to remove and collect large volumes of fluid. Tubing volume = approximately 25 mL from SENSAT.R.A.C.™ Pad to canister.

**Spinal Cord Injury:** In the event a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate in response to stimulation of the sympathetic nervous system), discontinue Negative Pressure Therapy to help minimize sensory stimulation.

**Bradycardia:** To minimize the risk of bradycardia, the dressing must not be placed in proximity to the vagus nerve.

**Enteric Fistula or Leak:** When treating an open abdomen where enteric fistulas are present, clinicians should consider the potential for abdominal contamination if effluent is not appropriately isolated or managed.


**Protect Periwound Skin:** Consider use of a skin preparation product to protect periwound skin. Do not allow foam to overlap onto intact skin. Protect fragile / friable periwound skin with additional drape, hydrocolloid or other transparent film.

- Multiple layers of the drape may decrease the moisture vapor transmission rate, which may increase the risk of maceration.
- If any signs of irritation or sensitivity to the drape, foam or SENSAT.R.A.C.™ Pad tubing appear, discontinue use and consult a physician.
- To avoid trauma to the periwound skin, do not pull or stretch the drape over the foam dressing during drape application.

If there are any questions regarding the proper placement or usage of the ABTHERA ADVANCE™ Open Abdomen Dressing, please contact your local KCI clinical representative.

### DRESSING APPLICATION

#### ABTHERA ADVANCE™ OPEN ABDOMEN DRESSING COMPONENTS

### WOUND PREPARATION


**WARNING: Review all ABTHERA ADVANCE™ Open Abdomen Dressing Safety Information before beginning wound preparation. Ensure adequate hemostasis has been achieved prior to dressing placement (refer to Bleeding section under WARNINGS).**

1. Sharp edges or bone fragments must be eliminated from wound area or covered (refer to **Bleeding** section under **WARNINGS**).
2. Irrigate abdominal wound and cleanse periwound skin as indicated.
3. Clean and dry periwound tissue; consider use of a skin preparation product to protect periwound skin. Do not allow foam to overlap onto intact skin. Protect fragile / friable periwound skin with additional drape, hydrocolloid or other transparent film.

### VISCERAL PROTECTIVE LAYER APPLICATION


The Visceral Protective Layer is fenestrated to allow for active fluid removal when negative pressure is applied and is designed to allow application of this layer directly over omentum or exposed internal organs.

**WARNING: The foam in the Visceral Protective Layer is encapsulated for patient safety. Protect vital structures with Visceral Protective Layer at all times during therapy. Never place exposed foam material directly in contact with exposed bowel, organs, blood vessels or nerves.**


1. Remove contents from inner pouch and unfold the Visceral Protective Layer in a sterile field. Either side of the Visceral Protective Layer may be placed on the omentum or viscera.
2. Gently place Visceral Protective Layer over the open abdominal cavity (**Fig. 1**).
3. Determine the orientation of the dressing for the specific application. If Visceral Protective Layer will be placed around tubes, drains or the falciform ligament, cut only between the foam extensions (**Fig. 2**). Do not cut near or through foam extensions. Orient the Visceral Protective Layer accordingly before cutting.
4. Size the Visceral Protective Layer by folding or cutting as described in the following sections.

#### Folding Visceral Protective Layer to Size

1. Hold dressing by the edge and slightly lift. Slowly lower dressing into the paracolic gutter, using the other hand to gently and evenly work the dressing down (**Fig. 3**). Fold any excess Visceral Protective Layer up and over onto itself.
2. Continue placing Visceral Protective Layer between abdominal wall and internal organs (**Fig. 4**) throughout the abdominal compartment. The goal is to ensure full coverage of all viscera.


#### Cutting Visceral Protective Layer to Size

1. Cut Visceral Protective Layer away from wound, through center of large foam squares using sterile scissors (**Fig 5A**). Do not cut through narrow connecting tabs between the large foam squares.
2. Pinch the remaining half of the foam square and its connecting tab and pull. The foam and tab will separate at the next square (**Fig. 5B**). This will ensure that edges of Visceral Protective Layer cover exposed foam edge (**Fig. 5C**) and foam cannot contact organs (see **WARNING** above).
3. Document number of foam extensions removed and that each piece has been properly disposed of away from wound cavity.

**CAUTION:** Do not tear the foam over the wound, as fragments may fall into the wound. Rub or trim foam away from wound, removing any fragments to ensure loose particles will not fall into or be left in the wound upon dressing removal.

## PERFORATED FOAM APPLICATION

The perforated foam (Fig. 6) provided with the ABTHERA ADVANCE™ Open Abdomen Dressing is intended to:

- Transfer negative pressure from the Negative Pressure Therapy Unit to the Visceral Protective Layer to promote active fluid removal.
- Provide medial tension upon foam collapse to help maintain fascial domain.


1. Tear or cut perforated foam to needed size as shown above (Fig. 7). The foam should fit directly over the Visceral Protective Layer and be in contact with wound edges. Do not allow foam to contact intact skin. One or both pieces of the provided perforated foam can be used, depending on the wound profile.

2. Gently place perforated foam into wound cavity over the Visceral Protective Layer (Fig. 8).

**NOTE:** Ensure foam-to-foam contact for even distribution of negative pressure.

**NOTE:** Always note the total number of pieces of foam used and document on the drape and in the patient's chart.

## DRAPE APPLICATION

1. Holding the drape, partially pull back one side of layer 1 to expose adhesive (Fig. 9). Be sure to hold layer 1 flap back, to prevent re-adherence to drape.

2. Place the drape adhesive-side down to cover foam and intact skin, ensuring drape covers at least an 8 - 10 cm border of intact periwound tissue (Fig. 10). Use any excess drape to seal difficult areas, if needed.

**NOTE:** To avoid trauma to the periwound skin, do not pull or stretch the drape over the foam dressing. Minimize wrinkles, as they may be a source of negative pressure leaks (refer to **PRECAUTIONS, Protect Periwound Skin** section).


3. Remove remaining tab 1 backing material and pat around drape to ensure an occlusive seal.

4. Remove green-striped stabilization layer 2 (Fig. 11).

5. Remove perforated blue handling tabs from drape (Fig. 12).

**NOTE:** When using multiple pieces of drape, ensure that the edges of the drape overlap in order to achieve a seal (Fig. 13).

## SENSAT.R.A.C.™ PAD APPLICATION


**NOTE:** Do not cut off the pad or insert the tubing into the foam dressing. This may occlude the tubing and cause the Negative Pressure Therapy Unit to alarm and could injure underlying viscera.

1. Choose pad application site. Give particular consideration to fluid flow and tubing position to allow for optimal flow and avoid placement over bony protuberances or within creases in the tissue.
2. Pinch drape and cut a 2.5 cm hole (not a slit) through the drape (Fig. 14). It is not necessary to cut into the foam.

**NOTE:** Cut a hole rather than a slit, as a slit may self-seal during therapy.

3. Apply pad, which has a central disc and a surrounding outer adhesive skirt.
  - Gently remove both backing layers 1 and 2 to expose adhesive (Fig. 15).
  - Place pad opening in central disc directly over hole in drape (Fig. 16).
  - Apply gentle pressure on the central disc and outer skirt to ensure complete adhesion of the pad.
4. Pull back on blue tab to remove pad stabilization layer (Fig. 17). Dressing application is complete.

## V.A.C.® NEGATIVE PRESSURE THERAPY APPLICATION

**NOTE:** Only for use with Negative Pressure Therapy provided by INFOV.A.C.™ and V.A.C.ULTA™ Negative Pressure Therapy Units. Refer to the therapy unit user manual for complete instructions for use.

**NOTE:** SENSAT.R.A.C.™ Pad tubing is not compatible with hospital vacuum systems.

**WARNING: Review all Negative Pressure Therapy System Safety Information before initiating therapy.**

1. Remove canister from packaging and insert into the therapy unit until it locks into place.

**NOTE:** Abdominal wounds often have copious drainage. Consider using the 1000 cc / mL canister. Ensure an adequate supply of canisters is readily available.

**CAUTION:** Consider the size and the weight of the patient, patient condition, wound type, monitoring capability and care setting when using the 1000 cc / mL canister.

**NOTE:** If the canister is not fully engaged, the therapy unit will alarm.

2. Connect SENSAT.R.A.C.™ Pad tubing to canister tubing and ensure clamp on each tube is open (Fig. 18). Position clamps away from patient.

3. Turn on power to the therapy unit and select 125 mmHg, continuous mode therapy setting for efficient fluid removal rates. Negative pressure therapy settings below 125 mmHg are not recommended.

**CAUTION:** Do not use intermittent therapy / DYNAMIC PRESSURE CONTROL™ Therapy with the ABTHERA ADVANCE™ Open Abdomen Dressing.

4. Initiate therapy. Assess dressing to ensure integrity of seal. The dressing should collapse and have a wrinkled appearance. There should be no hissing sounds. If there is any evidence of non-integrity, check drape and SENSAT.R.A.C.™ Pad seals, tubing connections, and canister insertion, and ensure clamps are open. Secure excess tubing to prevent inadvertent tension on tubing, which may disrupt the seal.

**Monitor Fluid Output** - The dressing is designed to efficiently remove fluid from the abdominal compartment and to evenly distribute negative pressure. When treating patients with the Negative Pressure Therapy Unit, the volume of exudate in the canister and tubing should be frequently examined

**Bleeding: Patients with abdominal wounds must be closely monitored for bleeding as these wounds may contain hidden blood vessels which may not be readily apparent. If sudden or increased bleeding is observed in the dressing, tubing or canister, immediately discontinue Negative Pressure Therapy, take appropriate measures to stop bleeding, and contact the physician. Negative Pressure Therapy is not designed to prevent, minimize or stop bleeding. (Refer to WARNINGS, Bleeding section).**

## ALARM RESOLUTIONS

All therapy unit alarms should be addressed in a timely manner. Refer to the therapy unit user manual for complete information on alarm resolutions.

In case of a leak alarm, patch leak source with additional drape to ensure integrity of seal.

**CAUTION:** Due to the highly exudative nature of abdominal wounds, Negative Pressure Therapy should be interrupted only for wound care and dressing change. Interruption of therapy can result in loss of seal integrity.

## DRESSING CHANGES

Dressing changes should occur every 24 to 72 hours, or more frequently based upon a continuing evaluation of wound condition and patient presentation. Consider more frequent dressing changes in the presence of infection or abdominal contamination.

Refer to **Application Setting** section under **WARNINGS**.

Whenever the dressing is changed, always replace all dressing components with components from an unopened sterile package.

## DRESSING REMOVAL

Remove and discard previous dressing per institution protocol. Completely inspect wound, including paracolic gutters, to ensure all pieces of dressing components have been removed. If intra-abdominal packing is present, packing material may be drier than anticipated. Evaluate packing material prior to removal and rehydrate if necessary to prevent adherence or damage to adjacent structures.

**WARNING: Refer to Dressing Removal section under WARNINGS.**

## EXPLANATION OF SYMBOLS USED

	Do not use if package is damaged or open		Manufacturer
	Single use only		Date of Manufacture
	Do not reutilize		Content information
	Consult instructions for use		
	Contains Phthalates		
	DEHP		
	Method of sterilization - Radiation		Catalog number
	Use by		Lot number
	Keep dry		

**Rx Only CAUTION:** Federal (US) law restricts this device to sale / rental by or on the order of a physician.

Always count and record number of foam pieces used in wound.

References available on request. Please contact KCI at 1-800-275-4524 (in the US).

Kaplan M. Managing the open abdomen. Ostomy Wound Management, 2004 Jan; 50(1A suppl); C2, 1-8

Kaplan M, Banwell P, Orgill DP, Ivatury RR, Demetriades D, Moore FA, Miller P, Nicholas J, Henry S. Guidelines for the Management of the Open Abdomen. WOUNDS. 2005 Oct; 17(Suppl 1); S1S24

Garner GB, Ware DN, Cocanour CS, Duke JH, McKinley BA, Kozar RA, Moore FA. Vacuum-assisted wound closure provides early fascial reapproximation in trauma patients with open abdomens. The American Journal of Surgery, 2001 Dec; 182(6); 630-8

Barker DE, Kaufman HJ; Vacuum Pack Technique of Temporary Abdominal Closure; A 7-Year Experience with 112 Patients. Presented at the 59th Annual Meeting of the American Association for the Surgery of Trauma. September 16-18, 1999. Boston Mass.

Brock WB, Barker DE; Temporary Closure of Open Abdominal Wounds; The Vacuum Pack. Presented at the 66th Annual Scientific Meeting of the Southeastern Congress, Lake Buena Vista, Florida. February 6-10, 1994

Sherck J, Seiver A; Covering the "Open Abdomen"; A Better Technique. Presented as a Poster at the 66th Annual Scientific Meeting and the Postgraduate Course Program. Southeastern Surgical Congress. Atlanta, Georgia. January 31-February 4, 1998.

## CONTACT INFORMATION

For questions regarding this product, maintenance, or additional information about KCI products and services, please contact KCI or a KCI authorized representative, or:

In the US call 1-800-275-4524 or visit [www.aclecity.com](http://www.aclecity.com) or [www.openabdomen.com](http://www.openabdomen.com). KCI USA, Inc., 12930 IH 10 West, San Antonio, TX 78249

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Garner GB, Ware DN, Cocanour CS, Duke JH, McKinley BA, Kozar RA, Moore FA. Vacuum-assisted wound closure provides early fascial reapproximation in trauma patients with open abdomens. The American Journal of Surgery, Dez. 2001; 182(6); 630-8

Barker DE, Kaufman HJ; Vacuum Pack Technique of Temporary Abdominal Closure; A 7-Year Experience with 112 Patients. Vorgestellt bei der Veranstaltung 59th Annual Meeting of the American Association for the Surgery of Trauma. 16.–18. September 1999. Boston Mass.

Brock WB, Barker DE; Temporary Closure of Open Abdominal Wounds; The Vacuum Pack. Vorgestellt bei der Veranstaltung 66th Annual Scientific Meeting of the Southeastern Congress, Lake Buena Vista, Florida. 6.–10. Februar, 1994

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Garner GB, Ware DN, Cocanour CS, Duke JH, McKinley BA, Kozar RA, Moore FA. Vacuum-assisted wound closure provides early fascial reapproximation in trauma patients with open abdomens. The American Journal of Surgery, 2001 Dec; 182(6); 630-8.

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Brock WB, Barker DE; Temporary Closure of Open Abdominal Wounds; The Vacuum Pack. Presented at the 66th Annual Scientific Meeting of the Southeastern Congress, Lake Buena Vista, Florida. February 6-10, 1994.

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Kaplan M, Barwell P, Orgill DP, Ivatury RR, Demetriades D, Moore FA, Miller P, Nicholas J, Henry S, Guidelines for the Management of the Open Abdomen. WOUNDS. 2005 Oct; 17(Suppl 1); S1524

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Barker DE, Kaufman HJ; Vacuum Pack Technique of Temporary Abdominal Closure: A 7-Year Experience with 112 Patients. Presented at the 59th Annual Meeting of the American Association for the Surgery of Trauma. September 16-18, 1999. Boston Mass.

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Garner GB, Ware DN, Cocanour CS, Duke JH, McKinley BA, Kozar RA, Moore FA. Vacuum-assisted wound closure provides early fascial reapproximation in trauma patients with open abdomens. The American Journal of Surgery, diciembre 2001; 182(6); 630-8

Barker DE, Kaufman HJ; Vacuum Pack Technique of Temporary Abdominal Closure; A 7-Year Experience with 112 Patients. Presentado en la 59.ª Annual Meeting of the American Association for the Surgery of Trauma. 16-18 septiembre, 1999. Boston Mass.

Brock WB, Barker DE; Temporary Closure of Open Abdominal Wounds; The Vacuum Pack. Presentado en el 66.º Annual Scientific Meeting of the Southeastern Congress, Lake Buena Vista, Florida. 6-10 febrero, 1994

Sherck J, Seiver A; Covering the «Open Abdomen»; A Better Technique. Presentado como póster en el 66.º Annual Scientific Meeting y el Postgraduate Course Program. Southeastern Surgical Congress. Atlanta, Georgia. 31 de enero al 4 de febrero, 1998.

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Copyright 2018 KCI Licensing, Inc. Alle rechten vorbehalten. 419020 Rev A 3/2018
Direitos autorais 2018 KCI Licensing, Inc. Todos os direitos reservados. 419020 Rev A 3 / 2018