Epidermis-Entnahmesystem Benutzerhandbuch / Epidermaal oogststelsysteem Gebruikershandleiding

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Suomi

Norsk

Epidermal Harvesting System
User Manual

IMPORTANT SAFETY INFORMATION ACCOMPANIES THIS DEVICE

Warnings, Precautions and other Safety Information are contained in this CelluTome™ Epidermal Harvesting System User Manual included with the control unit. Please consult this user manual before starting a procedure.

If there are questions, contact your local KCI representative.

Additional product information can be found at www.kci1.com (USA) or www.kci-medical.com (outside the USA).

As with all prescription medical devices, failure to follow warnings, precautions and instructions for use may lead to improper product performance and the potential for patient injury. For medical questions please consult a physician.

CAUTION: Federal law (US) restricts this device to sale or rental by or on the order of a physician.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

KCI HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE KCI PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY KCI SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL KCI BE LIABLE FOR ANY INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND KCI TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in KCI printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties except as set forth in the written limited warranty included with this product. Information in this publication may be subject to change at any time. Contact KCI for updates.

WARNINGS: IMPORTANT INFORMATION FOR USERS

In order for KCI products to perform properly, KCI recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with this manual and applicable product labeling.
- Do not modify the control unit or vacuum head. Do not connect the control unit or vacuum head to other devices that may be in use.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by KCI. For these authorized personnel, KCI will make available upon request circuit diagrams, component parts lists, etc. as required for repairs.
- Ensure the electrical installation of the room complies with the appropriate national electrical wiring standards. To avoid the risk of electrical shock, this product must be connected to a grounded power receptacle.
- Do not operate this product if it has a damaged power cord or plug. If these components are worn or damaged, contact KCI.
- Do not drop or insert any object into any opening or tubing of this product.
- Do not connect this product or its components to devices not recommended by KCI.
- Only use accessories and harvesters provided with this product.
- Keep this product away from heated surfaces.
- Although this product conforms to the intent of the standard IEC 60601-1-2 in relation to Electromagnetic Compatibility, electrical equipment, including cell phones or similar products, may produce interference. If interference is suspected, separate the equipment and contact KCI.
- Avoid spilling fluids on any part of this product.

WARNING: Fluids remaining on the electronic controls can cause corrosion that may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and staff. If spills do occur, unplug the unit immediately and clean with an absorbent cloth.

- Ensure there is no moisture in or near the power connection and power supply components before reconnecting power. If the product does not work properly, contact KCI.
- Refer to the **Care and Cleaning** section of this manual for information on infection control.

Notice

This product has been configured from the manufacturer to meet specific voltage requirements. Refer to the product information label for specific voltage.

TABLE OF CONTENTS

Important Safety Information Accompanies This Device	2
Disclaimer of Warranty and Limitation of Remedy.....	2
Warnings: Important Information For Users.....	3
Intended Use.....	5
Contraindications.....	5
Device Description	5
Warnings	5
Precautions.....	5
Control Unit User Interface Panel.....	6
Preoperative Setup.....	7
Vacuum Hose Connection.....	7
Power Up	7
Harvesting Instructions.....	8
Control Unit Troubleshooting Guide.....	10
Control Unit Stand Assembly Instructions.....	11
Tools Required:.....	11
Mounting the Control Unit onto the Stand.....	12
Care and Cleaning.....	13
Standard Precautions.....	13
Waste Disposal.....	13
Cleaning the Vacuum Head	13
Cleaning the Control Unit.....	13
Cleaning the Control Unit stand.....	13
Maintenance.....	13
Reordering Information.....	14
CelluTome™ Epidermal Harvesting System - Components and Accessories.....	14
CelluTome™ Epidermal Harvesting System - Disposables.....	14
Electromagnetic Compatibility.....	15
CelluTome™ Epidermal Harvesting System - Power Cord Information.....	19
Explanation of Symbols Used	20
Specifications.....	21
Customer Contact Information	22

INTENDED USE


The CelluTome™ Epidermal Harvesting System is intended to reproducibly cut a thin skin graft for autologous skin grafting.

CONTRAINDICATIONS

None


DEVICE DESCRIPTION

The **reusable** components of the CelluTome™ Epidermal Harvesting System:

- The control unit creates and regulates the vacuum and monitors the vacuum head.
- The vacuum head and tubing deliver the vacuum from the control unit and creates the warming required to raise the epidermal suction microdomes.

The **disposable** component of the CelluTome™ Epidermal Harvesting System:

The harvester holds the suction microdomes that are raised for procurement. A 3M™ Tegaderm™ Film (catalog number: 1624W) is placed onto the harvester's microdome grid prior to cutting the microdomes. After cutting the microdomes, the 3M™ Tegaderm™ Film with the micrografts is then applied to the recipient site.

Follow warnings, precautions and instructions listed in this document. Improper use could result in patient injury, bleeding, scarring, inadequate micrografts, prolonged procedure time, incomplete or ineffective procedure, the need to perform harvesting at a second donor site or patient discomfort.

WARNINGS

The harvester is a disposable component and is intended for single use only. Re-use of disposable components may result in contamination and infection.

PRECAUTIONS

Inspect the CelluTome™ Epidermal Harvesting System before use for signs of visible damage. Do not use if there are any indications of damage.

Harvest from areas of healthy intact skin (donor site).

Actuating the harvester handle prior to complete microdome formation may result in patient pain and / or bleeding.

Actuating the harvester handle prior to placement of the 3M™ Tegaderm™ Film may result in an inability to secure the micrografts for application.

Wipe the vacuum head and the control unit nest following institutional procedures used for cleaning and disinfection of other hard surface durable electronic medical equipment between each patient use. Inadequate cleaning may result in patient contamination.

Use this system only with the supplied components.


Discard all disposable items (harvester and 3M™ Tegaderm™ Film) in accordance with local medical waste disposal regulations.

Contact KCI for disposal of the control unit at the end of its expected service life. This product is designated for separate collection at an appropriate collection point. Do not dispose of this product in normal waste stream.

Handle the vacuum head with care. Do not drop the vacuum head, as doing so could damage the glass warming window. Do not suspend or hold the vacuum head by the tubing

Reversible maceration at the recipient site is a potential risk. Please use best clinical judgement for preventing maceration. Fenestration of the 3M™ Tegaderm™ Film with an 18g needle prior to application onto the recipient site may be appropriate.

CONTROL UNIT USER INTERFACE PANEL

Status Indicator:

- steady amber - unit plugged in
- steady green - unit turned on
- flashing green - unit running
- flashing amber - fault condition (see below)

Vacuum Fault: Indicates that sufficient vacuum to perform the procedure has not been achieved. In this case ensure that the harvester is securely strapped onto the patient's donor site, and the vacuum head is securely snapped onto the harvester.


Other Fault / Service Unit: Do not use the system. Contact your KCI representative.

NOTE: Refer to page 10 for control unit interface panel messages and troubleshooting.

PREOPERATIVE SETUP

VACUUM HOSE CONNECTION

1. Locate the tubing connections inside the cradle of the control unit.
2. Plug the vacuum head's heater tube (A) and the vacuum tube (B) into the control unit connections as shown below.

POWER UP

1. Plug control unit into a power outlet.

To isolate the control unit from supply mains, unplug the AC power from the wall outlet.


Position the unit so that the power cord connection is accessible at all times to allow for immediate disconnection from power source, if necessary.

2. Press On / Off switch on the outside of the control unit's lower back panel. The Status Indicator light will be amber.
3. Press the **Power** button on the user interface panel to turn the control unit on. The Status Indicator light will change from amber to green.

HARVESTING INSTRUCTIONS

The harvester is provided sterile.

1. Prepare donor site by shaving or clipping hair in a 10 cm square at the donor site if needed. Rigorously wipe area with alcohol.
2. Open the disposable harvester package.
3. With the blue handle oriented up, position the harvester on the patient's inner thigh, centered on the area of prepared skin, then secure with the integrated strap. Confirm complete contact with donor site by ensuring skin can be seen in all of the microdome grid holes in the harvester. Reposition if necessary.
4. Ensure that the vacuum head tubing is plugged into the control unit, the control unit is plugged into the wall and that the Status Indicator light is amber.
5. Press the **Power** button on the user interface panel to turn the control unit on. The Status Indicator light will change from amber to green.
6. Press the blue latches to detach the vacuum head and remove it from the control unit storage nest.
7. With the tubing facing up, place the vacuum head onto the harvester, and push until the vacuum head snap fits onto the harvester. Ensure that the vacuum head is securely latched to the harvester.
8. The timer will initially display 30, indicating that the timer will beep after 30 minutes (Default Alert Time).

NOTE: The control unit does not automatically shut off at 30 minutes. After 30 minutes, the timer will beep every minute until it reaches 75 minutes. After 75 minutes the control unit will automatically shut off.

9. Use the **Timer Up** and **Timer Down** buttons to adjust the Default Alert Time as desired.

NOTE: The Default Alert Time can only be adjusted in 5 minute increments (minimum of 5, maximum of 75).

10. Press **Start / Pause** button to begin the suction microdome raising process. The Status Indicator light will flash green. If needed, press the **Start / Pause** button again to momentarily pause the suction microdome raising process.
11. Once the timer beeps, check the microdomes to see if they are all ready to harvest - raised, clear, and shiny. (Fig. 1)
12. If the system was paused and microdomes are not ready to harvest, press the **Start / Pause** button to continue the process.
13. Once the microdomes are ready to harvest, press the **Power** button or the **Start / Pause** button to stop the microdome raising process.
14. Unlatch the vacuum head from the harvester. The harvester should remain on the patient's skin.
15. Wipe the underside of the vacuum head following institutional procedures, then place the vacuum head in the control unit nest and push to engage the latches.
16. Place 3M™ Tegaderm™ Film onto the harvester grid. Center the film starting from the middle then fan out to the edges to ensure complete coverage of the protruding microdomes.
17. Use fingers to firmly press the 3M™ Tegaderm™ Film against microdomes. Apply firm pressure to ensure that the 3M™ Tegaderm™ Film adheres to the microdomes.
18. Harvest the microdomes by raising the blue harvester handle fully (until an audible click is heard) and then lowering the handle to the original position.
19. Carefully peel back the 3M™ Tegaderm™ Film from one end of the harvester. Use care to ensure that the 3M™ Tegaderm™ Film does not get folded, torn or otherwise damaged when removing it from the harvester. The 3M™ Tegaderm™ Film with attached micrografts is now ready for application and should be applied within two minutes of procurement.
20. Apply the 3M™ Tegaderm™ Film with micrografts to the recipient site.
21. Apply additional pressure dressings / bandages as appropriate for the recipient site per your current standard of care.
22. Remove the harvester from the patient's donor site.
23. Apply a second 3M™ Tegaderm™ Film or other treatment / dressing to the donor site per your current standard of care.



Fig. 1

CONTROL UNIT TROUBLESHOOTING GUIDE

The following section describes states of operation for the CelluTome™ Control Unit. If fault conditions cannot be resolved, contact your KCI representative.

Off (no power)

Standby – power available, user interface panel is not active

- Timer display is blank
- Status Indicator light is solid amber
- Vacuum and heater off
- LEDs off in vacuum head

Power On – user interface panel is active after **Power** button is pressed

- Timer shows default time (30 min)
- Status Indicator light is solid green
- Vacuum and heater off
- LEDs off in vacuum head

Running – system functioning, after **Start / Pause** button is pressed

- Timer shows running time elapsed (counting up from zero)
- Status Indicator light is flashing green
- Vacuum and heater on
- LEDs on in vacuum head

Paused – after start of procedure

- Timer flashes actual running time elapsed (not counting)
- Status Indicator light is solid green
- Vacuum and heater off
- LEDs off in vacuum head

Complete – pre-set run time has elapsed

- Timer continues counting
- Status Indicator light is flashing green
- Vacuum and heater on
- Timer beeps starting at time set and every minute after
- LEDs on in vacuum head

Complete / Stopped – 75 minutes has elapsed (safety feature)

- Timer shows flashing 75
- Status Indicator light is steady green
- Vacuum and heater off
- LEDs on in vacuum head

Vacuum Fault: – vacuum seal has been compromised (or never made)


- Timer flashes running time elapsed (not counting)
- Status Indicator light flashes amber after ten seconds and the **vacuum icon** is illuminated
- Vacuum and heater off
- LEDs off in vacuum head

Other Fault / Service Unit: – heater broken, tubing unplugged, any problem other than vacuum

- Timer shows two underscores (_ _)
- Status Indicator light flashes amber and the **service icon** is illuminated
- Vacuum and heater off
- LEDs off in vacuum head


CONTROL UNIT STAND ASSEMBLY INSTRUCTIONS

TOOLS REQUIRED:

- 1/2" Box Wrench
1. Insert base pole with tapered bushing into top of base.
 2. Screw hex head bolt with washers into base pole as shown.
 3. Tighten securely with 1/2" box wrench.

MOUNTING THE CONTROL UNIT ONTO THE STAND

1. Position the control unit as shown.
2. Ensure the location pins are aligned with the location holes.
3. Push the control unit down until it fits securely onto the base pole and location pins and is flush with the mount plate.

CARE AND CLEANING

KCI recommends the following regarding cleaning and disinfecting therapy devices:

- To help reduce risk of infection and contact with blood and body fluids, use personal protective equipment (PPE) such as medical procedure gloves.
- Do not immerse or saturate the control unit with fluids to avoid damage to the electronics in the device

STANDARD PRECAUTIONS

The following are the KCI recommended cleaning and infection control procedures for the CelluTome™ Epidermal Harvesting System.

Always follow Standard Precautions.

Standard Precautions are designed to reduce the risk of transmission of microorganisms from both known and unknown sources of infection. These precautions can be applied to all patients, regardless of their diagnosis or presumed infection status, and should be used when contact is anticipated with blood and all body fluids. This also includes secretions and excretions (except sweat) regardless of whether blood is visible or not, non-intact skin (i.e., open wounds) and mucous membranes.

WASTE DISPOSAL

Discard all disposable items (harvester and 3M™ Tegaderm™ Film) in accordance with local medical disposal regulations.

CLEANING THE VACUUM HEAD

After each use, wipe down the vacuum head following institutional procedures used for cleaning and disinfection of other hard surface durable electronic medical equipment.

CLEANING THE CONTROL UNIT

As needed, wipe down the control unit following institutional procedures used for cleaning and disinfection of other hard surface durable electronic medical equipment.

Ensure that the control unit and its power supply are not connected to AC power when using cleaning fluids of any nature.

CLEANING THE CONTROL UNIT STAND

Wipe the stand with a damp towel and a mild cleaning solution once a week, or more often if required.

NOTE: *If IV solution spills on base, clean immediately.*

MAINTENANCE

There is no maintenance required on the part of the user.

REORDERING INFORMATION

CELLUTOME™ EPIDERMAL HARVESTING SYSTEM - COMPONENTS AND ACCESSORIES

Catalog Number	Description	Quantity
CT-KIT/US	CelluTome™ Control Unit and Vacuum Head with Carrying Case for US and Canada	Kit, One Each
CT-KIT/EU	CelluTome™ Control Unit and Vacuum Head with Carrying Case for EU	Kit, One Each
CT-KIT/GB	CelluTome™ Control Unit and Vacuum Head with Carrying Case for UK and Ireland	Kit, One Each
CT-KIT/IT	CelluTome™ Control Unit and Vacuum Head with Carrying Case for Italy	Kit, One Each
CT-KIT/DK	CelluTome™ Control Unit and Vacuum Head with Carrying Case for Denmark	Kit, One Each
CT-KIT/CH	CelluTome™ Control Unit and Vacuum Head with Carrying Case for Switzerland	Kit, One Each
CT-KIT/AU	CelluTome™ Control Unit and Vacuum Head with Carrying Case for AU and NZ	Kit, One Each
CT-KIT/ZA	CelluTome™ Control Unit and Vacuum Head with Carrying Case for South Africa	Kit, One Each
CT-CUS1	CelluTome™ Control Unit Stand	Each

CELLUTOME™ EPIDERMAL HARVESTING SYSTEM - DISPOSABLES

Catalog Number	Description	Quantity
CT-H25	CelluTome™ Harvester, 2.5 cm	Five per box
CT-H50	CelluTome™ Harvester, 5.0 cm	Five per box
CT-TFLM	CelluTome™ 3M™ Tegaderm™ Film	100 per box

ELECTROMAGNETIC COMPATIBILITY

Although this equipment conforms to the intent of the 2004/108/EC EMC Directive, all medical equipment has the potential to cause electromagnetic interference or be susceptible to electromagnetic interference. The following are guidance and manufacturer's declarations regarding EMC for the CelluTome™ Epidermal Harvesting System.

- The CelluTome™ Epidermal Harvesting System needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following pages.

WARNING: This equipment is intended for use by healthcare professionals only. As with all electrical medical equipment, this equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures such as re-orienting or relocating the CelluTome™ Epidermal Harvesting System or shielding the location.

- Portable and Mobile RF communications equipment, RFID readers, electronic article surveillance (anti-theft) equipment and metal detectors can affect the performance of the CelluTome™ Epidermal Harvesting System. Please use the guidelines and recommendations in the following pages.
- Other medical equipment or systems can produce electromagnetic emissions and therefore can interfere with the functionality of the CelluTome™ Epidermal Harvesting System. Care should be used when operating the CelluTome™ Epidermal Harvesting System adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the CelluTome™ Epidermal Harvesting System should initially be observed to verify normal operation in the configuration in which it will be used.
- The electrical cables, external power supplies and accessories listed or referenced in this manual have been shown to comply with the test requirements listed in the following tables. Care should be taken to use only manufacturer recommended cables, power supplies and accessories with the CelluTome™ Epidermal Harvesting System. If a third-party supplier offers cables, external power supplies and electrical accessories for use with the CelluTome™ Epidermal Harvesting System and they are not listed or referenced in this manual, it is the responsibility of that third-party supplier to determine compliance with the standards and tests in the following tables.
- The use of electrical cables and accessories other than those specified in this manual or referenced documents may result in increased electromagnetic emissions from the CelluTome™ Epidermal Harvesting System or decreased electromagnetic immunity of the CelluTome™ Epidermal Harvesting System.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The CelluTome™ Epidermal Harvesting System is intended for use in the electromagnetic environment specified below. The customer or the end user of the CelluTome™ Epidermal Harvesting System should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The CelluTome™ Epidermal Harvesting System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	The CelluTome™ Epidermal Harvesting System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CelluTome™ Epidermal Harvesting System is intended for use in the electromagnetic environment specified below. The customer or the end user of the CelluTome™ Epidermal Harvesting System should assure it is used only in such an environment.

IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input / output lines	±2 kV for power supply lines ±1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CelluTome™ Epidermal Harvesting System requires continued operation during power mains interruptions, it is recommended that the CelluTome™ Epidermal Harvesting System 'be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the A.C. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CelluTome™ Epidermal Harvesting System is intended for use in the electromagnetic environment specified below. The customer or the user of the CelluTome™ Epidermal Harvesting System should assure that it is used in such an environment.

IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the CelluTome™ Epidermal Harvesting System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation Distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CelluTome™ Epidermal Harvesting System is used exceeds the applicable RF compliance level above, the CelluTome™ Epidermal Harvesting System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CelluTome™ Epidermal Harvesting System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the CelluTome™ Epidermal Harvesting System

The CelluTome™ Epidermal Harvesting System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CelluTome™ Epidermal Harvesting System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CelluTome™ Epidermal Harvesting System as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter W	Separation Distance According to Frequency of Transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

CELLUTOME™ EPIDERMAL HARVESTING SYSTEM - POWER CORD INFORMATION

Power Cord	Description	Cord Specification	Max Length (inches)
350084	Cord, AC Power	3 x 18 AWG, SJT, 10A / 125V	78.74
350753	Cord, V.A.C. ULTA™ Power, UK-240V	H05VVF-3G, 10A / 250V	78.74
360074	Cord, V.A.C. VIA™ Power, UK-220V	H05VVF-3G, 10A / 250V	79.00
360076	Cord, V.A.C. VIA™ Power, AU / NZ-240V	H05VVF-3G, 10A / 250V	79.00
360080	Cord, V.A.C. VIA™ Power, IT-220V	H05VVF-3G, 10A / 250V	79.00
350758	Cord, V.A.C. ULTA™ Power, DK-220V	H05VVF-3G, 10A / 250V	78.74
360081	Cord, V.A.C. VIA™ Power, CH-220V	H05VVF-3G, 10A / 250V	79.00
360122	Cord, V.A.C. VIA™ Power South Africa	H05VVF-3G, 10A / 250V	79.00

WARNING: The use of electrical cables and accessories other than those specified in this manual or referenced documents may result in increased electromagnetic emissions from the CelluTome™ Epidermal Harvesting System or decreased electromagnetic immunity of the Cellutome™ Epidermal Harvesting System.

EXPLANATION OF SYMBOLS USED

Method of Sterilization -
Ethylene Oxide

Keep Dry

Do not use if package
is damaged or open

Manufacturer

Single Use Only

This product is designated for separate collection at
an appropriate collection point. Do not dispose of
as household waste.

Do Not Resterilize

Non-sterile

Rx Only

CAUTION: Federal (US) law restricts this device
to sale / rental by or on the order of a physician.

Consult Instructions
For Use

IP2

Protection against solid objects greater than 12.5 mm

Refer to User Manual

Type BF
Applied Part

Conforms with the Medical Device
Directive (93/42/EEC) and has been
subject to the conformity procedures
laid down in the council directive.

Temperature
Limitations

Fragile

ETL Listed, Conforms to AAMI ES60601-1 1st
edition, CSA C22.2#60601-1 3rd edition and
IEC 60601-1 3rd edition

Date of Manufacture

Content Information

Lot Number

REF

Catalog Number

Authorized Representative in
the European Community

Use By

SPECIFICATIONS

Specifications subject to change without notice.

Harvester

Dimensions:5.7"W x 3.7"H x 6"D (14.5 x 9.4 x 15.2 cm)

Weight:~0.2 lbs (~0.09 kg)

Vacuum Head

Dimensions:5.5"W x 3.9"H x 1.3"D (13.9 x 9.9 x 3.3 cm)

Weight:~0.6 lbs (~0.27 kg)

Control Unit

Dimensions:10.6"W x 9.3"H x 9.4"D (26.9 x 23.6 x 23.9 cm)

Weight:~7.3 lbs (~3.31 kg)

Control Unit Basket Accessory

Maximum weight rating:2 lbs (the approximate weight of 4 individual harvester packages)

Electrical

Internal Power Supply Input:100-240 VAC 3.2A 50 / 60 Hz

Internal Power Supply Output:24 VDC, 10.4 A

Power Inlet Fuse

- Type: 5x20mm Time Lag glass cartridge
- Current Rating: 5A
- Voltage Rating: 250V
- Operative Speed(s): 600ms ≤ t ≤ 10s @ 13.75A
- Breaking Capacity: 150A @ 250V

Storage / Transport Conditions

Temperature Range:0°F (-18°C) to 140°F (60°C)

Relative Humidity Range:15 - 95% non-condensing

Operating Conditions

Temperature Range:58°F (14°C) to 95°F (35°C)

Relative Humidity Range:20 - 80% non-condensing

Atmospheric Pressure:700 hPa to 1013 hPa

Expected Service Life:3 years

IEC Classification

- Medical Equipment
- Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- Type BF, Applied Part
- Class I
- IP2
- The harvester components of the CelluTome™ Epidermal Harvesting System are considered Applied Parts under IEC 60601-1 Third Edition.

CUSTOMER CONTACT INFORMATION

For questions regarding this product, supplies, maintenance, or additional information about KCI products and services, please contact KCI or a KCI authorized representative, or:

In the US call 1-800-275-4524 or visit www.kci1.com.

Outside the US visit www.kci-medical.com.

KCI USA, Inc.
12930 IH 10 West
San Antonio, Texas 78249 USA
www.kci1.com

KCI Medical Products (UK), Ltd.
11 Nimrod Way
Wimborne, Dorset
BH21 7SH
United Kingdom
www.kci-medical.com

Rx Only

3M™ and Tegaderm™ are trademarks of 3M™ Corporation. All other trademarks designated herein are proprietary to KCI Licensing, Inc. and its affiliates and / or licensors.

©2014 KCI Licensing, Inc. All rights reserved. 415307 Rev B 03/2014

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©2014 KCI Licensing, Inc. Med enerett. 415307 Rev B 03/2014