SNAP™ SecurRing™ HYDROCOLLOID

INSTRUCTIONS FOR USE

ENGLISH

For use exclusively with the SNAP™ Therapy System

PRODUCT DESCRIPTION

The SNAP™ SecurRing™ Hydrocolloid is a sterile, moldable hydrocolloid ring, designed to make dressing applications faster and easier. The Hydrocolloid can be molded into uneven skin surfaces or challenging body contours, filling these voids and reducing potential leak paths. It will also absorb moisture and help protect intact skin from excessive body fluids.

Not made with natural rubber latex.

INDICATIONS FOR USE

The SNAP™ SecurRing™ Hydrocolloid is indicated for use:

- To fill uneven skin surfaces, scars, skin wrinkles, and skin folds to help create a negative pressure seal.
- To mold around anatomical contours reducing potential leak paths
- To protect intact skin from excessive moisture


CONTRAINDICATIONS

- Do not use on patients who have a sensitivity or allergic reaction to hydrocolloids.

WARNINGS

- For external use only.
- Do not apply in an open wound or skin ulcer.
- Should irritation occur, discontinue use and consult a health care professional.
- Do not use on patients who are allergic to adhesives.
- For single patient use only. Do not reuse, reprocess, or resterilize.
- The SNAP™ SecurRing™ Hydrocolloid is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken or difficult to open.
- Discard product if mishandling has caused possible damage or contamination.
- Do not use if SNAP™ SecurRing™ Hydrocolloid is damaged or contaminated.
- During removal of the SNAP™ SecurRing™ Hydrocolloid, ensure that all pieces of hydrocolloid are removed from the wound area.
- The SNAP™ SecurRing™ Hydrocolloid should be stored at room temperature only and should not be exposed to excessive cold or heat.

APPLICATION INSTRUCTIONS

1. Make sure skin is clean, dry, and free from any solvent or greasy substances.
2. Remove the SNAP™ SecurRing™ Hydrocolloid from packaging (**Fig. 1**). Remove the release liner from both surfaces.
3. Gently mold, stretch, compress, or roll the SNAP™ SecurRing™ Hydrocolloid to any shape to fit around the outside of the patient's wound (**Fig. 2**). The SNAP™ SecurRing™ Hydrocolloid may be pulled apart and rejoined to better fit around the patient's wound. For smaller wounds, a portion of the SNAP™ SecurRing™ Hydrocolloid can be used.

4. Gently apply to intact skin around the outside of the wound with light pressure (**Fig. 3**). Ensure there is at least 1 cm of intact skin between the wound border and SNAP™ SecurRing™ Hydrocolloid.

NOTE: Do not place SNAP™ SecurRing™ Hydrocolloid directly into wound.

5. Follow the SNAP™ Therapy System instructions for use included with dressings and dressing kits to apply contact layer and dressing on wound.

6. Apply the SNAP™ Therapy Dressing over contact layer and SNAP™ SecurRing™ Hydrocolloid per normal procedures.

NOTE: Ensure dressing seals completely over SNAP™ SecurRing™ Hydrocolloid. Do not allow SNAP™ SecurRing™ Hydrocolloid to extend beyond the edges of the dressing.

SNAP™ SecurRing™ HYDROCOLLOID REMOVAL


1. Gently peel the SNAP™ Therapy Dressing from the body. The SNAP™ SecurRing™ Hydrocolloid should peel off with the dressing.
2. Gently roll and peel any remaining SNAP™ SecurRing™ Hydrocolloid to remove from the skin.
3. Ensure all pieces of SNAP™ SecurRing™ Hydrocolloid are removed from wound area.

DISPOSAL

Adhere to institution disposal protocol and requirements.

For detailed instructions and product information on the SNAP™ Therapy System which includes warnings and precautions, please refer to the Instructions for Use enclosed in the SNAP™ Therapy Dressing Kits.

EXPLANATION OF SYMBOLS USED

KCI 0473
KCI USA, Inc.
12930 IH 10 West
San Antonio, TX 78249 USA
1-800-275-4524
www.acelity.com

EC REP
KCI Medical Products (UK), Ltd.
11 Nimrod Way
Wimborne, Dorset
BH21 7SH
United Kingdom
www.kci-medical.com

KCI
An Acelity Company

All trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and/or licensors.
©2016 KCI Licensing, Inc. All rights reserved.
L22754 Rev D 6/2016

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Alle andre varemerker som nevnes her, tilhører KCI Licensing, Inc., deres datterselskaper og/eller lisensgivere.
©2016 KCI Licensing, Inc. Med enerett.