REVOLVE™

ADVANCED ADIPOSE SYSTEM

Table of Contents

| | |
|---|-------|
| Instructions for Use English | 2-4 |
| Instrucciones de uso Español | 5-7 |
| Mode d'emploi Français..... | 8-11 |
| Gebrauchsanweisung Deutsch | 12-15 |
| Istruzioni per l'uso Italiano..... | 16-19 |
| Instruções de utilização Português..... | 20-23 |
| Gebruiksaanwijzing Nederlands..... | 24-27 |
| Brugsanvisning Dansk | 28-30 |
| Instruktioner för användning Svenska | 31-33 |
| Bruksanvisning Norsk..... | 34-36 |
| Notkunarleiðbeiningar Íslenska | 37-40 |
| Instruções de uso Português do Brasil | 41-44 |

Device Description and Intended Use:

REVOLVE™ System is a sterile single-use disposable tissue canister used for harvesting, filtering, separating, concentrating, and transferring autologous tissue components for reintroduction to the same patient during a single surgical procedure for repair, reconstruction, or replacement of integumentary or musculoskeletal tissues. For additional volume of tissue, canister may be filled and emptied a second time during the same surgery on the same patient.

Contraindications:

Contraindications to autologous fat transfer include the presence of any disease processes that adversely affect wound healing, and poor overall health status of the individual.

Warnings:

1. Reuse of REVOLVE™ System may result in infection and the transmission of communicable diseases.
2. If sterile packaging is damaged do not use product.
3. This device will not, and is not intended to, produce significant weight reduction.
4. This device should be used with extreme caution in patients with chronic medical conditions such as diabetes, heart, lung, or circulatory system disease or obesity.
5. The volume of blood loss and endogenous body fluid loss may adversely affect intra and/or postoperative hemodynamic stability and patient safety. The capability of providing adequate, timely replacement is essential for patient safety.


Precautions:

1. Use of this device is limited to those physicians who, by means of formal professional training or sanctioned continuing medical education (including supervised operative experience), have attained proficiency in suction lipoplasty and tissue transfer.
2. Results of this procedure will vary depending upon patient age, surgical site, and experience of the physician.
3. Results of this procedure may or may not be permanent.
4. The amount of fat removed should be limited to that necessary to achieve a desired effect.

Adverse Effects:

Some common adverse effects associated with use of the REVOLVE™ System and/or autologous fat transfer procedures are asymmetry, over- and/or under-correction of the treatment site, tissue lumps, bleeding, scarring, fat necrosis, cyst formation, allergic reaction, and infection and inflammation of various levels. If an unanticipated event occurs, alteration of surgical plan may be necessary at the surgeon's discretion.

REVOLVE™ System Components


REVOLVE™ System Unit Parts

User Provided Components

1. Liposuction Tube

2. Connect to: Medical grade vacuum source with clamp.

3. Waste Canister

1/4" post-remove blue elbow
3/8" post-use blue elbow

4. Catheter-tip Syringe

or Luer Lock Syringe


5. Lactated Ringer's Solution

Step 1: Harvest

1. Leave REVOLVE™ System in tray.
2. Connect Vacuum Port to Waste Canister with connecting tube.
3. Connect Liposuction Tube to elbow (yellow) on Patient Port.
4. Hold device firmly to prevent tipping.
5. Harvest tissue into REVOLVE™ System using liposuction. Minimum harvest volume inside of mesh 100 ml.


Connect to medical grade vacuum source.

Step 2: Wash, Filter, and Concentrate

1. Disconnect/clamp Liposuction Tube from Patient Port.
2. Disconnect Vacuum at Vacuum Port OR turn off suction pump.
3. Open Vent Port. Leave open for the rest of the procedure.
4. Apply temperature strip to confirm Lactated Ringer's Solution is between 37-39°C. Connect only warm (37-39°C) Lactated Ringer's Solution to Patient Port using irrigation tubing. **MAINTAIN TEMPERATURE UP TO AND THROUGHOUT WASHING PROCESS.**
5. Fill Canister with approximately 300-350ml of warm (37-39°C) Lactated Ringer's Solution.
6. Rotate handle quickly for at least 15 seconds.
7. Reconnect Vacuum to Vacuum Port OR turn on suction pump, and remove fluid below mesh.
8. Repeat steps 5-7 two more times.
9. After third wash, maintain vacuum for 60 seconds to remove all fluids.


4. 37-39°C Lactated Ringer's Solution. Fluid must be warmed.

Temperature Strip

Note: Temperature strips will provide a visual indication of the Lactated Ringer's Solution temperature after it has been warmed to the appropriate temperature. Do not apply the temperature strip before warming the solution.

1. Remove liner from self-adhesive backing and ensure that the application surface of the Lactated Ringer's Solution is clean and dry before applying.
2. Green indicates the temperature. If green is not visible, the temperature will be midway between the illuminated blue and tan windows.
3. The temperature will only display if the Lactated Ringer's Solution is between 35°C-41°C.

Step 3: Transfer—IMMEDIATELY after Step 2

A. Catheter-tip Syringe

(Recommended method of extraction.)


1. **KEEP VENT PORT OPEN.**
2. Insert Catheter-tip Syringe into Patient Port.
3. Tilt Canister toward syringe, gently shake Canister to slide tissue toward syringe, then extract tissue.

Patient Port

If a syringe becomes clogged during extraction, push syringe plunger toward canister to release clog and then continue to extract.


4. Transfer to luer lock syringe via disposable adapter.

Warning: Do not overtighten luer connections

B. Luer Lock Syringe

1. Replace cap on Patient Port, open Vent Port.
2. **KEEP VENT PORT OPEN.**
3. Attach Luer Lock Syringe to Luer Port.
4. Tilt canister toward Syringe, gently shake Canister to slide tissue toward syringe, then extract tissue.

Luer Port

Guidelines for Tissue Temperature Management:

1. Use ONLY pre-warmed solutions for washing/filtering, warmed to temperature range of 37°C-39°C, per IFU. Cooler fluids cause the adipose to thicken and congeal. If room temperature fluids are used to wash, or the adipose is allowed to cool for more than a few minutes, the congealed adipose may clog the mesh filter. Utilizing solutions at temperatures below or above 37°C-39°C may lead to sub-optimal procedure outcomes.
2. Solutions can be pre-warmed to 37°C-39°C range in a fluid warmer or a warming chamber.
Alternately, a one liter bag can be pre-warmed in a microwave –
 - approximately 90 seconds using an 800 watt microwave;
 - approximately 72 seconds using a 1000 watt microwave;
 - approximately 60 seconds using a 1200 watt microwave;
 - approximately 48 seconds using a 1500 watt microwave.
 - Surgeons should confirm their own institutions' policies on use of microwaving prior to proceeding with these guidelines.
 - They should also assure that any microwaved Lactated Ringer's solution is in the 37°C-39°C range prior to inserting it in the REVOLVE™ System.
 - Once warmed, keep the solutions in a warmer to maintain 37°C -39°C range until filling canister. Ensure appropriate temperature prior to use.
3. **Begin the washing phase IMMEDIATELY after completion of the harvest phase. Do not let the tissue cool in the canister.**

Guidelines for Tissue Harvest:

1. Harvest at least 100 ml of adipose within the mesh filter. The maximum harvest is 350 ml of adipose within the mesh; do not overfill the canister, or the mesh filter may clog.
2. Harvest 20% - 30% more adipose than the desired or targeted amount. The washing phase reduces the volume of the adipose as harvested by washing out oils, interstitial fluids, trapped fluid pockets, and captures and retains collagen and other fibrous tissue.
3. For additional volume of tissue, canister may be filled and emptied a second time on same patient during same surgery.

Guidelines for Tissue Transfer:

1. Immediately extract the washed adipose into catheter tip syringes after washing is complete. If fat is thick or difficult to extract using catheter tip syringes, add a small amount of warm fluid to the tissue, stir, and then extract.
2. Use catheter tip syringes for tissue transfer from the device; Toomey syringes will not work.

Catheter tip syringes are legally, commercially available from all major syringe manufacturers. Common sizes are 30/35 ml, 50/60 ml, and 100/120 ml.

Definitions

CAUTION: FEDERAL (US) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

DO NOT USE IF PACKAGE IS DAMAGED OR OPEN

NON-PYROGENIC

CONTENTS

Dates provided on product labels are in the following format: 4 digit year, 2 digit month, and 2 digit day (YYYY-MM-DD).

This product, including its packaging and components, is not made with natural rubber latex.

US and international patents pending.

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© LifeCell Corporation October 2017
LifeCell Corporation
IPAX, Inc.
KCI Medical Products, Ltd.
KCI Medical Australia PTY, Ltd.