DE

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4.1.3 Split-skin application 11

4.1.4 Wound coverage 11

4.2 MatriDerm as a temporary adhesion barrier 11

4.3 MatriDerm for mucosal defects. 11

4.4 Other instructions for use. 12

5. Precautions 12

6. Contraindications 12

7. Side effects 12

8. Interactions 12

9. Shelf life and storage 12

10. Explanation of the symbols on the packaging 13

Collagen elastin template for tissue regeneration and use as a temporary adhesion barrier

1. Description

MatriDerm is a three-dimensional matrix composed of native structurally intact collagen fibrils and elastin to support tissue regeneration.

The collagen is obtained from bovine dermis and contains the dermal collagens I, III and V. The elastin is obtained from bovine nuchal ligament by hydrolysis.

MatriDerm serves as a scaffold for tissue regeneration and modulates the formation of scar tissue. Moreover, MatriDerm has excellent haemostatic properties and thus reduces the risk of subgraft haematoma, for example, between the matrix and the wound bed. In addition, MatriDerm supports orderly healing and acts as a temporary adhesion barrier.

The non-use of chemical cross-linking of the collagen results in a matrix with especially good biocompatibility. MatriDerm shows only limited mechanical stability when wetted, but still moulds very well to the wound bed.

2. Therapeutic indications

MatriDerm is used in combination with autologous split-thickness skin grafts for the reconstruction of deep dermal defects and full-thickness skin wounds in plastic reconstructive surgery and in the surgical treatment of burns, trauma and dermatological diseases. It is also used in the treatment of graft-requiring poorly healing wounds (e.g. chronic wounds). Additionally, MatriDerm can be used for the reconstruction of mucosal defects.

The aim of treatment is to construct a neodermis, in order to improve the quality of reconstituted tissue, reduce scarring, prevent wound contraction and restore functionality.

MatriDerm can also be used under intact skin for the temporary separation of tissues, thereby preventing adherence, especially of tendons and their surrounding connective tissue, after injuries or surgical procedures.

3. Mode of action

MatriDerm is an acellular tissue substitute, in which the native, structurally intact collagen acts as an

important constituent of the new extracellular matrix for the integration of cells and blood vessels. The low thickness of MatriDerm enables the graft to be initially supplied with nutrients by means of diffusion. Rapid vascularization then follows. As healing progresses, fibroblasts produce their own collagen matrix while MatriDerm is remodelled.

MatriDerm serves to promote an orderly tissue healing. After surgery, adhesences between various layers of tissue can develop due to unstructured scar formation during the healing process. Introduction of the acellular matrix provides the healing tissue with a scaffold. This prevents the formation of unstructured scar tissue that can lead to adherence. The movement of tissues against each other remains possible even after healing and consequential impairment caused by the adherence of tissues, such as restricted mobility and pain, is prevented.

The native structure of the collagen fibrils of MatriDerm has low antigenicity and outstanding biocompatibility.

In an animal study, MatriDerm gradually became fully degraded within four weeks of implantation. Experience from a human study on the treatment of punch biopsy wounds showed that MatriDerm was completely remodelled 6 weeks after implantation.

Investigations have shown that the elastin component of MatriDerm improves the stability and elasticity of the regenerating tissue. In an animal and human study, a reduction in wound contraction was also observed.

Clinical trials on the treatment of reconstructive wounds and burns showed that the elasticity of the regenerated skin after 3-4 months was significantly improved with the combined use of MatriDerm and split-thickness skin grafts than when the wounds were treated with split-thickness skin grafts alone. Long-term follow-up revealed that the regenerated skin of patients treated with MatriDerm still showed better elasticity 12 years after surgery.

4. Method of application

MatriDerm can be used as a dermal or mucosal substitute, or to avoid tissue adherence, especially

of the tendons. The method of use depends on the therapeutic indication and the recommendations are as follows:

4.1. MatriDerm as a dermal substitute

Depending on the treatment regime, MatriDerm can be applied in a one-step or two-step procedure. The one-step procedure can be used if the wound bed is capable of supplying nutrients to the overlying graft even through the 1 mm thick matrix by diffusion. The patient benefits in particular from the avoidance of a second operation.

For less well-vascularized wound beds or if ≥ 2 mm thick MatriDerm is necessary, then the two-step procedure should be used, because in these cases it cannot be safely assumed that the wound bed is capable of ensuring immediate and adequate nutrition of a split-thickness skin graft applied at the same time.

In both treatment regimes, the wound is first surgically debrided until a well-vascularized wound bed is present. Thorough haemostasis should be achieved before MatriDerm is used because an inadequate control of bleeding can result in the separation of MatriDerm; a residual haematoma jeopardizes integration of the overlying split-thickness skin.

4.1.1 One-step procedure

In the one-step procedure, MatriDerm is immediately covered with a split-thickness skin graft as described below.

This method gave good results, with early patient mobilization, in four independent clinical trials in burns patients and in one clinical trial with reconstructive surgery.

Note: Depending on the wound bed, a few days may be needed for the split-thickness skin graft to take.

Recommendation:

Dry application

Under sterile conditions, remove the appropriate piece of MatriDerm from the packaging and trim to roughly fit the skin defect. Now place the dry MatriDerm on the wound and initially press into the wound area using a laparotomy sponge. The MatriDerm starts to become moistened with wound secretion and adheres securely to the wound bed with no air bubbles. Trim the edges to a narrow circular overlap of approx. 2 mm.

Once the MatriDerm is correctly positioned in the wound area, use a syringe to carefully sprinkle saline onto the MatriDerm; alternatively, a very wet

laparotomy swab, soaked with saline, can be pressed onto MatriDerm and left in place for a few minutes. After this time, the MatriDerm is completely rehydrated and lies exactly in the wound bed corresponding to the prior placement.

Tiny white spots may remain in the material after rehydration in the wound area. However, this has proved to have no detrimental effect on graft take or on the post-operative results.

Alternatively:

Prior rehydration of MatriDerm

Rehydrate MatriDerm in plenty of physiological saline or Ringer solution before use.

Note: To avoid trapped pockets of air, MatriDerm should not be immersed in the solution or doused. Air pockets can hinder diffusion and thereby jeopardize the attached graft. Lay MatriDerm on the surface of the liquid so that it can become fully soaked within a few minutes.

The matrix is ready for use when the entire surface appearance has changed from "white" to "translucent".

Due to the mechanical properties of the rehydrated matrix, it is best to apply the matrix by hand. It is crucial that MatriDerm is in complete contact with the wound bed and adheres to it. Air bubbles between the wound bed and the matrix must therefore be carefully removed by smoothing them out to the margins of the matrix. If several matrices are required, position the individual pieces, so that they abut each other precisely or overlap slightly.

MatriDerm should be trimmed to match the wound edges.

Place the split-thickness skin graft in the wound area directly on top of MatriDerm. Sutures or staples can be used to provide additional anchoring of MatriDerm with the graft. If the matrix is fixed by suturing, take into account the limited tensile strength of the matrix when wet. Under certain circumstances, fixation can be achieved by tissue glue. If there is a delay in placement of the graft over the matrix, steps must be taken to ensure that MatriDerm does not dry out – suitable for this purpose is, for example, a sterile, non-adherent gauze soaked with physiological saline.

4.1.2 Two-step procedure

In the two-step procedure, MatriDerm is initially attached to the wound. As soon as good vascularization of the collagen-elastin template is achieved, final grafting with the split-thickness skin graft can take place in a second operation.

The skin can be grafted as soon as granulation tissue has grown through the matrix, typically after 5 to 12 days.

The granulation tissue that has extended beyond the matrix should be reduced to produce a suitable wound bed for the graft.

4.1.3 Split-skin application

The best results have been achieved with unmeshed skin grafts. Fenestration can help avoid seroma.

Experience in some cases suggests that good results can be achieved if the graft slightly overlaps the wound margins.

If a meshed skin graft is indicated, care must be taken to minimize the shear forces between the graft and MatriDerm. Unlike in the wound bed, adherence of the graft to MatriDerm is not promoted by the formation of fibrin.

4.1.4 Wound coverage

To ensure a moist wound environment, a non-adherent wound contact layer or a non-medicated tulle gras dressing is recommended as a wound dressing.

To date, it has been found practicable to use a conventional dressing (e.g. 3-4 sheets of gauze) or a negative pressure (vacuum) dressing over the top.

The use of a vacuum dressing on wounds treated with MatriDerm depends on the surgeon's preference in the particular patient. Good results have been achieved with this dressing technique in many cases.

Note: Do not use iodine-containing antiseptics in conjunction with MatriDerm. Caustic agents that alter proteins impair the function of collagen matrices and must not be used with MatriDerm.

The dressing technique applied must ensure good contact between the split-thickness skin graft, MatriDerm and the wound bed and prevent shear forces. Care must also be taken to ensure that the graft on MatriDerm cannot dry out, especially in the first few post-operative days.

Even with a correct dressing technique, isolated reports have appeared that the graft was still only weakly vascularized after 5 days, but took well a few days later.

4.2. MatriDerm as a temporary adhesion barrier

Recommendation:

Dry application

Under sterile conditions, remove the appropriate piece of MatriDerm from the packaging and trim to fit the target tissue. Be aware that in the dry state, the matrix exhibits wet adherence and can adhere to wet surfaces. When in the wetted adherent state, stability of the matrix to mechanical stress is only limited.

Place the dry MatriDerm carefully around or on the prepared tissue. When positioning the matrix in the operation site, it is recommended to allow the matrix to overlap the healthy tissue by a few millimeters.

In the case of cylindrical structures such as tendons or nerves, the matrix can be wrapped round once for fixing (twice the thickness in the overlap region).

Before normal closure of the wound, take care to ensure that the matrix is completely rehydrated. This can be achieved by the application of sterile physiological saline (alternatively, sterile Ringer-lactate solution) using a syringe. Depending on the wound, closure is achieved either by a direct method (e.g. with a suture, tissue adhesive or staples) or by covering MatriDerm with a split-thickness skin graft, as previously explained.

Transformation into endogenous tissue takes place within a few weeks.

4.3 MatriDerm for mucosal defects

MatriDerm can be used both in open as well as in closed, healing mucosal wounds.

Under sterile conditions, the wound is first surgically debrided until a well-vascularized wound bed is present. Thorough haemostasis should be achieved before MatriDerm is used because an inadequate control of bleeding can result in the separation of MatriDerm.

Remove the appropriate piece of MatriDerm from the packaging. Be aware that in the dry state, the matrix exhibits wet adherence and can adhere to wet surfaces. When in the wetted adherent state, stability of the matrix to mechanical stress is only limited.

Trim MatriDerm to fit the target tissue and place the dry MatriDerm carefully on the prepared tissue. When positioning the matrix in the operation site, it is recommended to allow the matrix to overlap the healthy tissue by a few millimeters. Take care to ensure that the matrix is completely rehydrated. This can be achieved by the application of e.g. sterile physiological saline.

Take care to remove all air bubbles from the matrix. If the matrix is fixed by suturing then take into account the limited tensile strength of the matrix when wet. Under certain circumstances, fixation can be achieved

by tissue glue. Common postoperative treatment measures and medication are usually sufficient.

4.4 Other instructions for use

MatriDerm is to be used only by surgeons.

When debriding a wound, care must be taken to completely remove scabs, necrosis and scar tissue.

MatriDerm is not indicated for use on devitalised tissue. However, a few cases of a successful grafting of MatriDerm onto poorly perfused wound beds with simultaneous covering with a split-thickness skin graft have been reported.

No noticeable reduction has been observed in the take rate of the graft onto MatriDerm compared with direct coverage of the wound base with a split-thickness skin graft. In one study with burn wounds, the growth rate of the graft in the one-step procedure was 83.4% on MatriDerm versus 82.5% without MatriDerm.

There are only a few studies in patients with chemical or electrical burns and/or radiation damage. The surgeon's decision to use MatriDerm for such wounds is based on the assessment of the wound and its suitability for surgical excision, the probable vitality of the wound bed and the balance of possible benefits for the patient compared to the risks.

5. Precautions

Because patients with radiotherapy and/or treatment with cytotoxics (e.g. chemotherapy) may suffer from wound healing disorders, special caution is also needed in these cases if MatriDerm is used.

Patients with untreated metabolic disorders (e.g. diabetes, osteomalacia and thyroid disease), long-term corticosteroid therapy, autoimmune conditions and heavy smoking habits may suffer from delayed wound healing. Therefore, in these patients, MatriDerm should only be used as a temporary adhesion barrier with special caution.

WARNING!

MatriDerm must not be resterilized!

MatriDerm must not be exposed to temperatures above 40°C, as otherwise the native collagen is irreversibly denatured.

6. Contraindications

Treatment with MatriDerm is contraindicated in patients with known hypersensitivity to bovine collagen or elastin.

Since a worsening of existing infections has been observed with other products, the use of MatriDerm in contaminated areas requires a particularly careful benefit-risk assessment by the treating surgeon.

The product must not be used in infected wounds!

No clinical data are available concerning the use of MatriDerm in pregnancy and during breastfeeding. Hence, the surgeon's decision to use MatriDerm must be based on consideration of all the risks in the particular case.

7. Side effects

Despite the low antigenic potential of collagen described in the literature, an intolerance reaction cannot be entirely excluded.

No cases of intolerance to MatriDerm have been reported to date, either in clinical trials with MatriDerm or during the continuous post-marketing surveillance.

General risks may include seromas, hematomas and hyperpigmentation. The patient should be made aware of these risks and others associated with this kind of surgery.

8. Interactions

Simultaneous use of collagenases, such as those used for enzymatic wound cleansing, can accelerate the degradation of MatriDerm.

Contact of the matrix with iodine-containing antiseptics or caustic agents that alter proteins is not indicated.

9. Shelf life and storage

The shelf life of MatriDerm (unopened in original packaging) is 60 months and it should not be used after the stated expiry date. The expiry date is imprinted on the packaging.

MatriDerm should be kept in a dry place at room temperature (do not store above 25° C).

MatriDerm is sterile-packed. To ensure sterility, intactness of the packaging should be checked. MatriDerm from opened sterile packs should be discarded. Do not resterilize MatriDerm, because this can affect the physicochemical properties of the matrix.

Keep medical devices out of the reach of children!

10. Explanation of the symbols printed on the package

Use by year + month

Upper temperature limit

Keep dry

Sterilized by using irradiation, sterile if unit package is undamaged and unopened

Do not re-use

Do not re-sterilize

Do not use if packaging is damaged

CE-sign with identification number of the Notified Body. The product complies with the essential requirements of Council Directive on Medical Devices 93/42/EEC.

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