plasma care®

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Hyggen (2020) Skin disinfection test according to VAH Method 13 on natural skin flora and supplemented with E. coli

Scheper et al. (2021) Cold plasma therapy with the handheld device plasma care® improves the tendency to heal in problem wounds - 10 case studies from diabetological practice

Yuta Terabe et. al. (2021) Using cold plasma to treat chronic foot ulcer infection

Yuta Terabe et. al. (2021) Treating hard-to-heal skin and nail onychomycosis of diabetic foot with plasma therapy

Brüning et. al. (2021) Using cold atmospheric plasma to treat hard to heal wounds. A case study with 10 Patients (with 19 wounds) treated in a dermatological outpatient clinic

Hämmerle [et.al.](#) 2021; Positive effects of cold atmospheric plasma on pH in wounds: a pilot study

Dejonckheere [et.al.](#) 2022. Non-Invasive Physical Plasma for Preventing Radiation Dermatitis in Breast Cancer: A First-In-Human Feasibility Stud

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Dejonckheere [et.al.](#) 2024. Non-invasive physical plasma for preventing radiation dermatitis in breast cancer: Results from an inpatient-randomised double-blind placebo-controlled trial

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² Isbary, G. et al. Atmospheric Plasma Devices for Medical Issues. Expert Rev. Med. Devices 2013, 10 (3), 367–377. <https://doi.org/10.1586/erd.13.4>.

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⁴ Isbary, G. et al. First Prospective Randomized Controlled Trial to Decrease Bacterial Load Using Cold Atmospheric Argon Plasma on Chronic Wounds in Patients. Br. J. Dermatol. 2010, 163 (1), 78–82. <https://doi.org/10.1111/j.1365-2133.2010.09744.x>.

⁵ Isbary, G. et al. Successful and Safe Use of 2 Min Cold Atmospheric Argon Plasma in Chronic Wounds: Results of a Randomized Controlled Trial. Br. J. Dermatol. 2012, 167 (2), 404–410. <https://doi.org/10.1111/j.1365-2133.2012.10923.x>.

⁶ Isbary, G. et al. Cold Atmospheric Argon Plasma Treatment May Accelerate Wound Healing in Chronic Wounds: Results of an Open Retrospective Randomized Controlled Study in Vivo. Clin. Plasma Med. 2013, 1 (2), 25–30. <https://doi.org/10.1016/j.cpm.2013.06.001>.

⁷ Heinlin, J. et al. Randomized Placebo- Controlled Human Pilot Study of Cold Atmospheric Argon Plasma on Skin Graft Donor Sites. Wound Repair Regen. Off. Publ. Wound Heal. Soc. Eur. Tissue Repair Soc. 2013, 21 (6), 800–807. <https://doi.org/10.1111/wrr.12078>.

⁸ Stratmann, B. et al. Effect of Cold Atmospheric Plasma Therapy vs Standard Therapy Placebo on Wound Healing in Patients With Diabetic Foot Ulcers: A Randomized Clinical Trial. JAMA Netw. Open 2020, 3 (7), e2010411. <https://doi.org/10.1001/jamanetworkopen.2020.10411>.

⁹ Jensen, J.-O. et al. The Repetitive Application of Cold Atmospheric Plasma (CAP) Improves Microcirculation Parameters in Chronic Wounds. Microvasc. Res. 2021, 138, 104220. <https://doi.org/10.1016/j.mvr.2021.104220>.

¹⁰ Arndt, S. et al. Effects of Cold Atmospheric Plasma (CAP) on β -Defensins, Inflammatory Cytokines, and Apoptosis-Related Molecules in Keratinocytes in Vitro and in Vivo. *PLoS One* 2015, 10 (3), e0120041. <https://doi.org/10.1371/journal.pone.0120041>.

¹¹ Arndt, S. et al. Cold Atmospheric Plasma (CAP) Activates Angiogenesis-Related Molecules in Skin Keratinocytes, Fibroblasts and Endothelial Cells and Improves Wound Angiogenesis in an Autocrine and Paracrine Mode. *J. Dermatol. Sci.* 2018, 89 (2), 181–190. <https://doi.org/10.1016/j.jdermsci.2017.11.008>.

¹² Arndt, S. et al. Cold Atmospheric Plasma (CAP) Changes Gene Expression of Key Molecules of the Wound Healing Machinery and Improves Wound Healing in Vitro and in Vivo. *PLoS One* 2013, 8 (11), e79325. <https://doi.org/10.1371/journal.pone.0079325>.

¹³ Hasse, S. et al. Induction of Proliferation of Basal Epidermal Keratinocytes by Cold Atmospheric-Pressure Plasma. *Clin. Exp. Dermatol.* 2016, 41 (2), 202–209. <https://doi.org/10.1111/ced.12735>.

¹⁴ Nicol, M. J. et al. Antibacterial Effects of Low-

Temperature Plasma Generated by Atmospheric-Pressure Plasma Jet Are Mediated by Reactive Oxygen Species. *Sci. Rep.* 2020, 10 (1), 3066. <https://doi.org/10.1038/s41598-020-59652-6>.

¹⁵ Hoon Park, J. et al. Comparative Study for the Inactivation of Multidrug Resistance Bacteria Using Dielectric Barrier Discharge and Nano-Second Pulsed Plasma. *Sci. Rep.* 2015, 5 (1), 13849. <https://doi.org/10.1038/srep13849>.

¹⁶ Zimmermann, J. et al. Test for Bacterial Resistance Build-up against Plasma Treatment. *New J. Phys.* 2012, 14 (7), 073037. <https://doi.org/10.1088/1367-2630/14/7/073037>.

¹⁷ Daeschlein, G. et al. Skin and Wound Decontamination of Multidrug-Resistant Bacteria by Cold Atmospheric Plasma Coagulation. *J. Dtsch. Dermatol. Ges. J. Ger. Soc. Dermatol. JDDG* 2015, 13 (2), 143–150. <https://doi.org/10.1111/ddg.12559>.

¹⁸ Bourke, P. et al. Microbiological Interactions with Cold Plasma. *J. Appl. Microbiol.* 2017, 123 (2), 308–324. <https://doi.org/10.1111/jam.13429>.

¹⁹ Weiss, Daeschlein, Kramer et al. Virucide Properties of Cold Atmospheric Plasma for Future Clinical Applications, *Journal of Medical Virology* 2017, 89, 952-959.

²⁰ Hämmerle et al. *Journal of Wound Care.* 2023, 32, 9, 530-536.

²¹ v. Brunn, Max-von-Pettenkofer Institut, unpublished data, showing that in first results plasma care can inactivate corona virus in solution

²² Lee, et al. Fast and easy disinfection of coronavirus-1 contaminated face masks using ozone gas produced by a dielectric barrier discharge plasma generator. *Medrxiv*, May 2020.

²³ Dargaville, T. et al. (2012). Sensors and Imaging for Wound Healing: A review. *Biosensors & bioelectronics.* 41. 10.1016/j.bios.2012.09.029.

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info@gdmedical.nl

www.gdmedical.nl

Jens Kirsch | Chief Executive Officer (CEO)

PD Dr. rer. nat. habil. med. Julia Zimmermann | Chief Medical Officer (CMO)

Prof. Dr. Dr. h.c. Gregor Morfill | Scientific Advisory Board

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terraplasma medical GmbH