(en)

## Peel and Place Dressing Kit

### Safety Information and Application Instructions

**Only for use with Solventum™ V.A.C.® Therapy Systems**

**Rx Only**

S solventum Internal Production Team	Solventum Production Artist: Brian Pedersen bpetersen@solventum.com	Division: MedSurg	This artwork has been produced by the Solventum Internal Production Team	Artwork Stage: 1	Date: 22-May-2025
File Name: 78-8145-4520-4b_IFU Solventum PandP Application - MDR.indb Structure Spec # N/A	Arwork ID: 0266840.1 Material Spec #: N/A Requestor: Anna Cantu PID #: 25702_g1				

FOR PRINTER USE ONLY: This template artwork was produced using a bar width adjustment of -0.20". If this does not meet the needs of the printing process or equipment to be used, the printer is responsible to recreate the file using the appropriate bar width adjustment while meeting Solventum quality requirements. Labels must be printed on a white background and must be printed in black ink. It is the responsibility of the printer to ensure that the label is legible and readable when applied to the final product. Should the supplied materials be found to be incorrect please contact a Solventum representative immediately. The liability of Solventum shall be limited to correcting its own product. Solventum will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost.

PMS Black

NON PRINT DATA

DIELINE

# PROOF REQUIRED

## **Table of Contents**

Solventum™ V.A.C.® Therapy System Safety Information.....	5
Indications for Use.....	5
Intended Users.....	5
Intended Patient Population.....	5
Contraindications.....	6
Warnings.....	6
Additional Warnings for the Solventum™ V.A.C.® Peel and Place Dressing.....	9
Precautions.....	10
Considerations for Transitioning Solventum™ V.A.C.® Therapy Into Home Care.....	12
Solventum™ V.A.C.® Peel and Place Dressing Kit Application Instructions.....	13
Component Identification.....	13
Solventum™ V.A.C.® Peel and Place Dressing Change Interval.....	14
Wound Preparation.....	14
Solventum™ V.A.C.® Peel and Place Dressing Removal.....	14
Solventum™ V.A.C.® Peel and Place Dressing Specific Information.....	15
Solventum™ V.A.C.® Peel and Place Dressing Application Instructions.....	16
Solventum™ SensaT.R.A.C.™ Pad Application.....	19
Solventum™ V.A.C.® Therapy Application with the Solventum™ V.A.C.® Peel and Place Dressing.....	20
Dressing Specifications.....	21
Symbols Glossary.....	22

## **Solventum™ V.A.C.® Therapy System Safety Information**

Disposable components of the Solventum™ V.A.C.® Therapy System are to be used as indicated on the associated product labeling. Solventum™ V.A.C.® Canisters are packaged sterile or fluid path sterile and are not made with natural rubber latex. Always use disposable components from sterile packages that have not been opened or damaged. All disposable components of the V.A.C.® Therapy System are for single use only. Re-use of disposable components may result in infection or impaired healing. To help ensure safe and effective use, the Solventum™ V.A.C.® Peel and Place Dressing Kit is to be used only with Solventum™ V.A.C.® Therapy Units.

The decision to use clean versus sterile / aseptic technique is dependent upon wound pathophysiology, physician / clinician preference and institutional protocol.

**Important:** As with any prescription medical device, failure to consult a physician and carefully read and follow all therapy unit and dressing instructions and safety information prior to each use may lead to improper product performance and the potential for serious or fatal injury. Do not adjust therapy unit settings or perform therapy application without directions from / or supervision by the treating physician.

### **Indications for Use**

The V.A.C.® Peel and Place Dressing Kit is intended to be applied in an acute, extended, or home care setting where product application is conducted by or under the supervision of a qualified healthcare professional.

When used on open wounds, the V.A.C.® Peel and Place Dressing Kit with the V.A.C.® Therapy System is intended to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. Open wound types are: acute and subacute dehisced wounds, traumatic wounds, partial-thickness burns, flaps, grafts and chronic diabetic ulcers, pressure ulcers and venous insufficiency ulcers.

### **Intended Users**

The V.A.C.® Peel and Place Dressing Kit is intended to be used by healthcare professionals in the acute, extended and home care settings. Dressing changes in the home care setting should be performed under the supervision of a healthcare professional.

### **Intended Patient Population**

The V.A.C.® Peel and Place Dressing Kit may be used on the general patient population as deemed appropriate by the prescribing clinician to promote wound healing in patients with wounds identified in the Indications for Use section.

## Contraindications

- Do not place dressings of the V.A.C.<sup>®</sup> Therapy System directly in contact with exposed blood vessels, anastomotic sites, organs or nerves.

**Note:** Refer to Warnings section for additional information concerning Bleeding.

- Solventum™ V.A.C.<sup>®</sup> Therapy is contraindicated for patients with:
  - Malignancy in the wound
  - Untreated osteomyelitis

**Note:** Refer to Warnings section for Osteomyelitis information.

- Non-enteric and unexplored fistulas
- Necrotic tissue with eschar present

**Note:** After debridement of necrotic tissue and complete removal of eschar, V.A.C.<sup>®</sup> Therapy may be used.

## Warnings

Please report a serious incident in relation to the device to Solventum and the local competent authority (EU) or the regulatory authority.

**Bleeding:** With or without using V.A.C.<sup>®</sup> Therapy, certain patients are at high risk of bleeding complications. The following types of patients are at increased risk of bleeding, which, if uncontrolled, could be potentially fatal:

- Patients who have weakened or friable blood vessels or organs in or around the wound as a result of, but not limited to:
  - Suturing of the blood vessel (native anastomoses or grafts) / organ
  - Infection
  - Trauma
  - Radiation
- Patients without adequate wound hemostasis
- Patients who have been administered anticoagulants or platelet aggregation inhibitors
- Patients who do not have adequate tissue coverage over vascular structures

If V.A.C.<sup>®</sup> Therapy is prescribed for patients who have an increased risk of bleeding complications, they should be treated and monitored in a care setting deemed appropriate by the treating physician.

If active bleeding develops suddenly or in large amounts during V.A.C.<sup>®</sup> Therapy, or if frank (bright red) blood is seen in the tubing or in the canister, immediately stop V.A.C.<sup>®</sup> Therapy, leave dressing in place, take measures to stop the bleeding and seek immediate medical assistance. The V.A.C.<sup>®</sup> Therapy Units and dressings should not be used to prevent, minimize or stop vascular bleeding.

- **Protect Vessels and Organs:** All exposed or superficial vessels and organs in or around the wound must be completely covered and protected prior to the administration of V.A.C.<sup>®</sup> Therapy. Always ensure that Solventum™ V.A.C.<sup>®</sup> Dressings do not come in direct contact with vessels or organs. Use of a thick layer of natural tissue should provide the most effective protection. If a thick layer of natural tissue is not available or is not surgically possible, multiple layers of meshed non-adherent material or bio-engineered tissue may be considered as an alternative, if deemed by the treating physician to provide a complete protective barrier. If using non-adherent materials, ensure they are secured in a manner that will maintain their protective position throughout therapy.

- Consideration should also be given to the negative pressure wound therapy setting and therapy mode used when initiating therapy. Caution should be taken when treating large wounds that may contain hidden vessels which may not be readily apparent. The patient should be closely monitored for bleeding in a care setting deemed appropriate by the treating physician.
- **Infected Blood Vessels:** Infection may erode blood vessels and weaken the vascular wall which may increase susceptibility to vessel damage through abrasion or manipulation. **Infected blood vessels are at risk of complications, including bleeding, which, if uncontrolled, could be potentially fatal. Extreme caution should be used when V.A.C.<sup>®</sup> Therapy is applied in close proximity to infected or potentially infected blood vessels.** (Refer to Protect Vessels and Organs section). The patient should be closely monitored for bleeding in a care setting deemed appropriate by the treating physician.
  - **Hemostasis, Anticoagulants, and Platelet Aggregation Inhibitors:** Patients without adequate wound hemostasis have an increased risk of bleeding, which, if uncontrolled, could be potentially fatal. These patients should be treated and monitored in a care setting deemed appropriate by the treating physician.  
Caution should be used in treating patients on doses of anticoagulants or platelet aggregation inhibitors thought to increase their risk for bleeding (relative to the type and complexity of the wound). Consideration should be given to the negative pressure setting and therapy mode used when initiating therapy.
  - **Hemostatic Agents Applied at the Wound Site:** Non-sutured hemostatic agents (for example, bone wax, absorbable gelatin sponge or spray wound sealant) may, if disrupted, increase the risk of bleeding, which, if uncontrolled, could be potentially fatal. Protect against dislodging such agents. Consideration should be given to the negative pressure setting and therapy mode used when initiating therapy.
  - **Sharp Edges:** Bone fragments or sharp edges could puncture protective barriers, vessels or organs, causing injury. Any injury could cause bleeding, which, if uncontrolled, could be potentially fatal. Beware of possible shifting in the relative position of tissues, vessels or organs within the wound that might increase the possibility of contact with sharp edges. Sharp edges or bone fragments must be covered or eliminated from the wound area to prevent them from puncturing blood vessels or organs before the application of V.A.C.<sup>®</sup> Therapy. Where possible, completely smooth and cover any residual edges to decrease the risk of serious or fatal injury, should shifting of structures occur. Use caution when removing dressing kit components from the wound so that wound tissue is not damaged by unprotected sharp edges.
  - **1000 mL Canister:** DO NOT USE the 1000 mL canister on patients with a high risk of bleeding or on patients unable to tolerate a large loss of fluid volume, including children and the elderly. Consider the size and weight of the patient, patient condition, wound type, monitoring capability and care setting when using this canister. This canister is recommended for acute care (hospital) use only.
    - **Solventum™ V.A.C.<sup>®</sup> Rx4 Therapy Unit:** The V.A.C.<sup>®</sup> Rx4 Therapy Unit provides four independent therapy channels that may accommodate either 500 or 1000 mL canisters.  
**When using multiple channels on multiple wounds, DO NOT USE the 1000 mL canister for patients at high risk for excessive fluid loss.**

**Infected Wounds:** Infected wounds should be monitored closely and may require more frequent dressing changes than non-infected wounds, dependent upon factors such as wound conditions and treatment goals. Refer to Dressing Change Interval section for details regarding dressing change frequency. As with any wound treatment, clinicians and patients / caregivers should frequently monitor the patient's wound, periwound tissue and exudate for signs of infection, worsening infection or other complications. Some signs of infection are fever, tenderness, redness, swelling, itching, rash, increased warmth in the wound or periwound area, purulent discharge or strong odor. Infection can be serious, and can lead to complications such as pain, discomfort, fever, gangrene, toxic shock, septic shock and / or fatal injury. Some signs or complications of systemic infection are nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucus membranes, disorientation, high fever, refractory and / or orthostatic hypotension or erythroderma (a sunburn-like rash). **If there are any signs of the onset of systemic infection or advancing infection at the wound site, contact the treating physician immediately to determine if V.A.C.<sup>®</sup> Therapy should be discontinued.** For wound infections relating to blood vessels, please also refer to the section titled **Infected Blood Vessels**.

**Osteomyelitis:** V.A.C.<sup>®</sup> Therapy should NOT be initiated on a wound with untreated osteomyelitis as wound healing may not progress. Consideration should be given to thorough debridement of all necrotic, non-viable tissue, including infected bone (if necessary) and appropriate antibiotic therapy.

**Circulatory Compromise:** V.A.C.<sup>®</sup> Therapy should not be initiated on a wound with untreated, underlying circulatory compromise (e.g., acutely ischemic wounds). Consideration should be given to both treatment of the underlying cause(s) and therapy settings before using the V.A.C.<sup>®</sup> Therapy System.

**Dressing Placement:** Always use V.A.C.<sup>®</sup> Peel and Place Dressings from sterile packages that have not been opened or damaged or are expired as this may lead to infection. Do not force dressings into any area of the wound, as this may damage tissue, alter the delivery of negative pressure or hinder exudate and dressing removal.

**V.A.C.<sup>®</sup> Peel and Place Dressings are radiolucent, not detectable on X-Ray.**

**Dressing Removal:** Gently remove the dressing from the wound. Regardless of treatment modality, disruption of the new granulation tissue during any dressing change may result in bleeding at the wound site. Minor bleeding may be observed and considered expected. However, patients with increased risk of bleeding, as described in the **Bleeding** section, have a potential for more serious bleeding from the wound site. **If significant bleeding develops, immediately discontinue the use of the V.A.C.<sup>®</sup> Therapy System, take measures to stop the bleeding, and do not remove the foam dressing until the treating clinician is consulted. Do not resume the use of the V.A.C.<sup>®</sup> Therapy System until adequate hemostasis has been achieved and the patient is not at risk for continued bleeding.**

**Keep V.A.C.<sup>®</sup> Therapy On:** Never leave a V.A.C.<sup>®</sup> Peel and Place Dressing in place without active V.A.C.<sup>®</sup> Therapy for more than two hours. Inactive therapy could lead to maceration, infection or impaired healing. If therapy is off for more than two hours, remove the old dressing and irrigate the wound. Either apply a new V.A.C.<sup>®</sup> Peel and Place Dressing from an unopened sterile package and restart V.A.C.<sup>®</sup> Therapy, or apply an alternative dressing at the direction of the treating physician.

**Acrylic Adhesive and Silicone Layer:** The V.A.C.<sup>®</sup> Peel and Place Dressing has an acrylic adhesive coating and a silicone layer, which may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives or silicone. If a patient has a known allergy or hypersensitivity to such materials, do not use the V.A.C.<sup>®</sup> Therapy System. If any signs of allergic reaction or hypersensitivity develop, such as redness, swelling, rash, urticaria or significant pruritus, discontinue use and consult a physician immediately. If bronchospasm or more serious signs of allergic reaction appear, seek immediate medical assistance.

**Defibrillation:** Remove the V.A.C.<sup>®</sup> Peel and Place Dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and / or patient resuscitation.

**Magnetic Resonance Imaging (MRI) - V.A.C.<sup>®</sup> Therapy Units:** The V.A.C.<sup>®</sup> Therapy Units are MR Unsafe. Do not take the therapy unit into the MR environment.

**Magnetic Resonance Imaging (MRI) - V.A.C.<sup>®</sup> Peel and Place Dressing:** V.A.C.<sup>®</sup> Peel and Place Dressings are MR Safe and can remain on the patient while in an MR environment. Never leave a V.A.C.<sup>®</sup> Peel and Place Dressing in place without active V.A.C.<sup>®</sup> Therapy for more than two hours (refer to Keep V.A.C.<sup>®</sup> Therapy On section). Inactive therapy could lead to maceration, infection or impaired healing.

**Hyperbaric Oxygen Therapy (HBO):** Do not take the V.A.C.<sup>®</sup> Therapy Unit into a hyperbaric oxygen chamber. The V.A.C.<sup>®</sup> Therapy Unit is not designed for this environment and **should be considered a fire hazard**. After disconnecting the V.A.C.<sup>®</sup> Therapy Unit, either (i) replace the V.A.C.<sup>®</sup> Peel and Place Dressing with another HBO compatible material during the hyperbaric treatment or (ii) cover the unclamped end of the V.A.C.<sup>®</sup> Tubing with dry gauze. For HBO therapy, the V.A.C.<sup>®</sup> Peel and Place Tubing must not be clamped. Never leave a V.A.C.<sup>®</sup> Peel and Place Dressing in place without active V.A.C.<sup>®</sup> Therapy for more than two hours (refer to Keep V.A.C.<sup>®</sup> Therapy On section). Inactive therapy could lead to maceration, infection or impaired healing.

**Respiratory Compromise:** To reduce the risk of respiratory compromise (asphyxiation or choking on small parts), consideration should be given to the patient's mental capacity, hearing function and visual acuity. Check that the tubing is positioned so as not to cause an obstruction or injury, e.g. asphyxiation. Do not wrap the tubing around the patient's neck and ensure all removable parts are secured in their appropriate working positions. Do not modify any of the medical devices or accessories in any way or use unspecified parts. Choking may result if patient swallows a small part that has become detached from the device or accessory.

## **Additional Warnings for the Solventum<sup>™</sup> V.A.C.<sup>®</sup> Peel and Place Dressing**

**Wound Fillers:** DO NOT use any wound fillers with the V.A.C.<sup>®</sup> Peel and Place Dressing. If a wound filler is required, select an alternative dressing to use to avoid maceration, infection, or impaired healing.

**Cutting the Foam:** Do not cut the foam section of the V.A.C.<sup>®</sup> Peel and Place Dressing to avoid maceration, infection or impaired healing.

**Wound Depth:** Do not use the V.A.C.<sup>®</sup> Peel and Place Small Dressing on wounds with a greater depth than 2 cm. Do not use the V.A.C.<sup>®</sup> Peel and Place Medium Dressing on wounds with a greater depth than 4 cm. Do not use the V.A.C.<sup>®</sup> Peel and Place Large Dressing on wounds with a greater depth than 6 cm. If the wound is deeper than 6 cm, use an alternative V.A.C.<sup>®</sup> Dressing to avoid maceration, infection or impaired healing.

**Tunnels and Undermining:** Do not place the V.A.C.<sup>®</sup> Peel and Place Dressing into a wound with any tunnels or with undermining greater than 2 cm. If clinician considers wound to be inappropriate for the V.A.C.<sup>®</sup> Peel and Place Dressing, please use an alternative V.A.C.<sup>®</sup> Dressing to avoid serious injury.

**Negative Pressure Settings and Continuous Therapy Only:** The V.A.C.<sup>®</sup> Peel and Place Dressing must only be used at pressure settings between 75 and 150 mmHg in Continuous mode.

**Instillation / Solventum<sup>™</sup> Veraflo<sup>™</sup> Therapy:** DO NOT use with Veraflo Therapy (Instillation) provided by the Solventum<sup>™</sup> V.A.C.<sup>®</sup> Ultra Therapy Unit. Instillation into the wound dressed with the V.A.C.<sup>®</sup> Peel and Place Dressing may result in fluid leaks which may cause maceration.

**Open Abdomen:** DO NOT USE the V.A.C.<sup>®</sup> Peel and Place Dressing over the open abdomen or with Solventum<sup>™</sup> AbThera<sup>™</sup> Open Abdomen Negative Pressure Therapy. Use over the open abdomen may result in an inability to maintain a negative pressure seal which may cause infection, maceration or impaired healing.

**Use over Umbilicus:** If the V.A.C.<sup>®</sup> Peel and Place Dressing covers the umbilicus, the umbilicus must first be filled with an antimicrobial petroleum gauze prior to dressing application to avoid the risk of infection.

## Precautions

**Inherent Surgical Risk:** Medical interventions may include inherent risks associated with the operative procedure, concomitant therapies and/or co-morbidities.

**Standard Precautions:** To reduce the risk of transmission of bloodborne pathogens, apply standard precautions for infection control with all patients, per institutional protocol, regardless of their diagnosis or presumed infection status. In addition to gloves, use gown and goggles if exposure to body fluids is likely.

**Continuous versus Intermittent / Solventum<sup>™</sup> Dynamic Pressure Control<sup>™</sup> Therapy:**

Continuous rather than intermittent / Dynamic Pressure Control Therapy is required.

**Patient Size and Weight:** The size and weight of the patient should be considered when prescribing V.A.C.<sup>®</sup> Therapy. Infants, children, certain small adults and elderly patients should be closely monitored for fluid loss and dehydration. Also, patients with highly exuding wounds or large wounds in relation to the patient size and weight should be closely monitored, as these patients have a risk of excessive fluid loss and dehydration. When monitoring fluid output, consider the volume of fluid in both the tubing and canister.

**Spinal Cord Injury:** In the event a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate in response to stimulation of the sympathetic nervous system), discontinue V.A.C.<sup>®</sup> Therapy to help minimize sensory stimulation and seek immediate medical assistance.

**Bradycardia:** To minimize the risk of bradycardia, V.A.C.<sup>®</sup> Therapy must not be placed in proximity to the vagus nerve.

**Enteric Fistulas:** Wounds with enteric fistulas require special precautions to optimize V.A.C.<sup>®</sup> Therapy. V.A.C.<sup>®</sup> Therapy is not recommended if enteric fistula effluent management or containment is the sole goal of therapy.

**Protect Periwound Skin:** The V.A.C.<sup>®</sup> Peel and Place Dressing incorporates a hybrid acrylic and silicone drape so use of a skin preparation product is not necessary.

- Multiple layers of drape may decrease the moisture vapor transmission rate leading to moisture accumulation over intact skin, which may increase the risk of maceration or infection.
- If any signs of irritation or sensitivity to the drape, foam, or tubing assembly appear, discontinue use and consult a physician.
- To avoid trauma to the periwound skin, do not pull or stretch the dressing during application. Excessive force on intact skin may lead to tissue damage.
- Extra caution should be used for patients with neuropath etiologies or circulatory compromise.

**Circumferential Dressing Application:** Avoid use of circumferential dressings except in the presence of anasarca or excessively weeping extremities, where a circumferential drape technique may be necessary to establish and maintain a seal. Extreme care should be taken not to stretch or pull the drape when securing it, but let it attach loosely and stabilize the edges with an elastic wrap, if necessary. When using circumferential drape applications, it is crucial to systematically and recurrently palpate distal pulses and assess distal circulatory status. If circulatory compromise is suspected, discontinue therapy, remove dressing, and contact the treating physician / clinician.

**Pressure Points:** Periodically assess and monitor the location of tubing connectors, caps, clamps or other rigid components to ensure they do not create inadvertent pressure points in relation to patient position as it may lead to serious tissue damage or impaired healing.

**V.A.C.<sup>®</sup> Therapy Unit Pressure Excursions:** In rare instances, tubing blockages with the V.A.C.<sup>®</sup> Therapy Unit may result in brief vacuum excursions to more than 250 mmHg negative pressure. Resolve alarm conditions immediately. Refer to the Therapy Unit User Guide or Manual or contact your Solvентum representative for additional information.

**In addition to these general warnings and precautions for V.A.C.<sup>®</sup> Therapy, additional warnings and precautions apply to certain specialty V.A.C.<sup>®</sup> Dressings and V.A.C.<sup>®</sup> Therapy Units. Please refer to the specific product instructions for use and labeling prior to application.**

## **Considerations for Transitioning Solventum™ V.A.C.® Therapy Into Home Care**

**Warning:** Patients with an increased risk of bleeding complications should be treated and monitored in a care setting deemed appropriate by the treating physician.

In addition to the contraindications, warnings and precautions for use of V.A.C.® Therapy, consider the following before prescribing V.A.C.® Therapy for use in the home care setting.

- **The Patient's Situation:**

- Clinical condition (adequate hemostasis and a low risk of active and / or large amounts of bleeding at the wound site)
- Home environment (patient or family member / caregiver able to read and understand safety labeling, able to respond to alarms, able to follow instructions for use)

- **The Patient's Wound:**

- Must be assessed for exposed vessels, anastomotic sites, organs, and nerves. All exposed or superficial vessels and organs in or around the wound must be completely covered and protected prior to the administration of V.A.C.® Therapy (refer to Protect Vessels and Organs in the Warnings section).

- **The V.A.C.® Canister Size:**

- The 1000 mL canister is NOT intended for use in the home.

- **Labeling:**


- The prescribing physician and health care clinician should be familiar with the V.A.C.® Therapy instructional materials that accompany the therapy unit into the home.
- The prescribing physician and / or healthcare clinician should carefully review these materials with the patient and patient's caregiver.
- Solventum offers in-service and training programs for use of V.A.C.® Therapy. Contact your local Solventum representative. In the U.S., call 1-855-423-6725 for scheduling.

If there are any questions regarding the proper placement or usage of V.A.C.® Therapy, please contact your local Solventum representative. If package is open, damaged or expired, do not use and contact your local Solventum representative.

Please report a serious incident in relation to the device to Solventum and the local competent authority (EU) or the regulatory authority. For additional and most current information, please see Solventum's website at [eIFU.Solventum.com](http://eIFU.Solventum.com).

## **Solventum™ V.A.C.® Peel and Place Dressing Kit Application Instructions**

### **Component Identification**

Solventum™ V.A.C.® Peel and Place Dressing

Solventum™ V.A.C.® Ruler

Solventum™  
SensaT.R.A.C.™ Pad

Disposable components of the Solventum™ V.A.C.® Therapy System including the Solventum™ V.A.C.® Peel and Place Dressing, Solventum™ SensaT.R.A.C.™ Pad and Solventum™ V.A.C.® Ruler are packaged sterile and are not made with natural rubber latex. Solventum™ V.A.C.® Canisters are packaged sterile or fluid path sterile and are not made with natural rubber latex. All disposable components of the V.A.C.® Therapy System are for single use only. To help ensure safe and effective use, the V.A.C.® Peel and Place Dressing is to be used only with Solventum™ V.A.C.® Therapy Units.

The decision to use clean versus sterile/aseptic technique is dependent upon wound pathophysiology, physician/clinician preference and institutional protocol.

Always consult a physician and review and follow Solventum™ V.A.C.® Therapy Safety Information and V.A.C.® Therapy Unit instructions prior to use.

### **Device Compatibility**

The V.A.C.® Peel and Place Dressing is compatible with the following Solventum therapy systems:

- Solventum™ V.A.C.® Ultia Negative Pressure Wound Therapy System
- Solventum™ V.A.C.® Rx4 Negative Pressure Wound Therapy System
- Solventum™ ActiV.A.C.™ Negative Pressure Wound Therapy System
- Solventum™ V.A.C.® Simplicity Negative Pressure Wound Therapy System
- Solventum™ Prevena™ Plus Negative Pressure Wound Therapy System

The V.A.C.® Ultia and V.A.C.® Rx4 Therapy Systems are integrated negative pressure wound management systems for use in acute care settings and other professional healthcare environments where product use is conducted by or under the supervision of a qualified healthcare professional.

The ActiV.A.C., V.A.C.® Simplicity and Prevena Plus 125 Negative Pressure Wound Therapy Systems are integrated negative pressure wound management systems for use in acute, extended and home care settings.

## **Solvotentum™ V.A.C.® Peel and Place Dressing Change Interval**

Wounds being treated with the V.A.C.® Therapy System should be monitored on a regular basis. In a monitored, non-infected wound, **V.A.C.® Peel and Place Dressings may be left in place for up to 7 days with the frequency adjusted by the clinician as appropriate.** Infected wounds must be monitored often and very closely. For these wounds, dressings may need to be changed more often; the dressing change intervals should be based on a continuing evaluation of wound condition and the patient's clinical presentation, rather than a fixed schedule.

## **Wound Preparation**

**Warning:** Review all V.A.C.® Therapy System Safety Information before beginning Wound Preparation.

1. Remove and discard previous dressing per institution protocol or local, regional and national laws and provisions. Thoroughly inspect wound to ensure all pieces of dressing components have been removed.

## **Solvotentum™ V.A.C.® Peel and Place Dressing Removal**

2. Gently remove an existing V.A.C.® Peel and Place Dressing according to the following procedure:
  - a. Raise the tubing connectors above the level of the therapy unit.
  - b. Close clamp on the dressing tubing to avoid contaminants entering the tubing which could lead to wound infection.
  - c. Disconnect canister tubing from dressing tubing.
  - d. Allow the therapy unit to pull the exudate in the canister tube into the canister, then close the clamp on the canister tubing.
  - e. Press THERAPY ON/OFF to deactivate the therapy unit. Wait 15-30 seconds to allow foam to decompress.
  - f. Gently remove the V.A.C.® Peel and Place Dressing from the wound - gently stretch the drape horizontally to release adhesive from the skin. Do not peel vertically to avoid tissue damage.

Warning: Refer to Dressing Removal section under Warnings.

- g. Discard disposables according to institutional or hospital regulations or state, local, regional and national laws and provisions.

**Note:** If dressing adheres to the wound, consider introducing sterile water or normal saline into the dressing, waiting 15-30 minutes, then gently remove the dressing from the wound.

**If the patient complains of discomfort during the dressing change, consider premedication as prescribed by the treating physician / clinician.**

3. Debride all necrotic, non-viable tissue, including bone, eschar or hardened slough, as prescribed by physician / clinician.
4. Perform thorough wound and periwound area cleaning per physician / clinician order or institution protocol prior to each dressing application.
5. Ensure adequate hemostasis has been achieved (refer to Warnings, Bleeding section, Hemostasis, Anticoagulants and Platelet Aggregation Inhibitors).

6. Prior to dressing placement, protect vessels and organs (refer to Warnings, Bleeding section, Protect Vessels and Organs).
7. Sharp edges or bone fragments must be eliminated from wound area or covered (refer to Warnings, Bleeding section, Sharp Edges).



### **Solventum™ V.A.C.® Peel and Place Dressing Specific Information**

- The V.A.C.® Peel and Place Dressing has a protective interface layer so the foam area can be placed over both the wound and surrounding intact skin. The interface layer allows the dressing to be used for up to 7 days.
- The drape border can be lifted and repositioned once during initial application to simplify dressing placement.
- Refer to the **V.A.C.® Peel and Place Dressing Application Instructions** section for appropriate **wound types and depths**. **An alternative V.A.C.® Dressing may be required for some wounds.**
- Refer to the **Solventum™ V.A.C.® Therapy Application with Solventum™ V.A.C.® Peel and Place Dressing** section for **recommended pressure settings**. Negative pressure must be sufficient to stretch the thin foam into the wound bed to avoid maceration, infection or impaired healing.
- The V.A.C.® Peel and Place Dressing is intended to be used **without sizing the drape border** for easy application, however, **trimming may be appropriate for some applications**, at the discretion of the clinician. **Leave adequate drape border for sealing. Failure to resolve leaks may cause infection, maceration or impaired healing.**

**Warning: Do not cut the foam section of the V.A.C.® Peel and Place Dressing to avoid maceration, infection or impaired healing.**

## Solventum™ V.A.C.® Peel and Place Dressing Application Instructions

1. Assess wound for any **undermining** or **tunneling**.


**Warning:** To avoid serious injury, do not place the V.A.C.® Peel and Place Dressing into a wound with **any tunnels or with undermining greater than 2 cm**. If clinician considers wound inappropriate for the V.A.C.® Peel and Place Dressing, please **use an alternative V.A.C.® Dressing**.


2. Assess wound dimensions and depth.

- **First**, choose dressing size based on **complete coverage** of the wound bed where NPWT is desired.


**Note:** The V.A.C.<sup>®</sup> Peel and Place Dressing foam should **extend beyond the edge of the wound** to ensure full coverage after V.A.C.<sup>®</sup> Therapy is applied.

- **Second**, check **maximum wound depth** recommended for the chosen V.A.C.<sup>®</sup> Peel and Place Dressing size. **If the wound is deeper than 6 cm, use an alternative V.A.C.<sup>®</sup> Dressing to avoid maceration, infection or impaired healing.**


3. Remove the release liners to expose adhesive.

4. Place the dressing onto the wound with the pre-cut hole facing up and positioned for the desired location of the SensaT.R.A.C. Pad. The hole does not need to be at the center of the wound.


**Note:** **Do not use wound fillers.** Use of wound fillers may lead to infection, maceration or impaired healing.

**Note:** To avoid trauma to the periwound skin, do not pull or stretch the dressing during application. Excessive force on intact skin may lead to tissue damage.


5. **Check for complete coverage of the wound where NPWT is desired.** The foam edges will be pulled inward when negative pressure is applied.

**Note:** If the dressing exposes wound edges after initial application, lift and readjust the dressing. **Do not lift and reapply during therapy. Lifting and reapplying dressing may result in disruption of wound area or infection.**

6. **Assess dressing placement.** If needed, the dressing **may be lifted**, repositioned and reapplied once for a better fit.

**Note:** Superficial or retention sutures should be covered with a single layer of non-adherent material placed between the sutures and the adhesive drape.

7. **Smooth** any wrinkles or creases in the dressing to prevent leaks. Failure to resolve leaks may cause infection, maceration or impaired healing.

**Note:** Drape can be peeled back and reapplied during initial placement to address affected part of the dressing.


8. While holding down the edge of the dressing, remove the perforated handling bars and pat down to ensure an occlusive seal.


## Solventum™ SensaT.R.A.C.™ Pad Application

**Note:** Do not cut off the pad or insert the tubing directly into the foam dressing. This may occlude the tubing and cause the V.A.C.® Therapy Unit to alarm.


1. The V.A.C.® Peel and Place Dressing comes pre-cut with a hole for the SensaT.R.A.C. Pad.

**Note:** Give particular consideration to tubing positioning to allow for optimal flow and avoid placement over bony prominences or within creases in the tissue to avoid serious tissue damage, or impaired healing.

2. Apply pad, which has a central disc and a surrounding outer adhesive skirt.
  - a. Remove both backing Layers 1 and 2 to expose adhesive.
  - b. Place pad opening in central disc directly over hole in drape.
  - c. Apply gentle pressure on the central disc and outer skirt to ensure complete adhesion of the pad.
  - d. Pull back on blue tab to remove pad stabilization layer.


## Solventum™ V.A.C.® Therapy Application with the Solventum™ V.A.C.® Peel and Place Dressing

**Warning:** Review all V.A.C.® Therapy System Safety Information before initiating V.A.C.® Therapy.

1. Remove canister from packaging and insert into the V.A.C.® Therapy Unit until it locks into place.

**Note:** If the canister is not fully engaged, the V.A.C.® Therapy Unit will alert / alarm.

2. Connect SensaT.R.A.C. Pad tubing to canister tubing and ensure clamp on each tube is open. Position clamps away from patient to avoid creating pressure points.
3. Turn on power to the V.A.C.® Therapy Unit and select the prescribed therapy setting.

**Warning:** The V.A.C.® Peel and Place Dressing **must only be used at negative pressure settings between 75 and 150 mmHg in Continuous mode to avoid maceration, infection or impaired healing.**

4. Initiate V.A.C.® Therapy. Assess dressing to ensure seal integrity. The dressing should be collapsed and have a wrinkled appearance. There should be no hissing sounds. For the ActiV.A.C., V.A.C.® Ultra and V.A.C.® Rx4 Therapy Systems, use the Solventum™ Seal Check™ Feature to verify that the rate of air leakage is below the alarm threshold. If there is any evidence of non-integrity, check SensaT.R.A.C. Pad and drape seals, tubing connections and canister insertion, and ensure clamps are open.

**Note:** Use caution when attaching device components as skin caught between connecting parts may lead to tissue damage.

If there is a leak, pat down the dressing or lift and reapply to ensure that there are no creases in the drape. Failure to resolve leaks may cause infection, maceration or impaired healing.

**Note:** Do not lift and reapply drape edges during therapy (only during initial placement).

5. Note the date that the dressing was applied and document on the supplied label on the V.A.C.® Ruler and in the patient's chart. The label can be peeled off the V.A.C.® Ruler and should be placed in an area that can be seen by the next treating clinician (placed around the SensaT.R.A.C. Pad tubing, in the patient's chart, etc.).

6. Secure excess tubing to prevent interference with patient mobility. Tubing may present a tripping hazard. Ensure that all tubing is out of areas where people may walk to reduce the risk of serious injury.

**Note:** Refer to unit specific user / clinician guide or manual and / or quick reference guide for information regarding alerts / alarms.

**Note:** If a leak source is identified, patch with additional drape to ensure seal integrity. Failure to resolve leaks may cause infection, maceration or impaired healing.

**Note:** If the wound is over a bony prominence or in areas where weight bearing may exert additional pressure or stress to the underlying tissues, a pressure redistribution (pressure relief) surface or device should be used to optimize patient offloading. Refer to Precautions, Pressure Point section.

**Warning:** Never leave a V.A.C.<sup>®</sup> Peel and Place Dressing in place without active V.A.C.<sup>®</sup> Therapy for more than two hours. Inactive therapy could lead to maceration, infection or impaired healing. If therapy is off for more than two hours, remove the old dressing and irrigate the wound. Either apply a new V.A.C.<sup>®</sup> Peel and Place Dressing from an unopened sterile package and restart V.A.C.<sup>®</sup> Therapy, or apply an alternative dressing at the direction of the treating physician / clinician.

## Dressing Specifications

Product	Overall Dressing Size	Foam Size	Maximum Wound Depth
EZ5SML - Solventum™ V.A.C. <sup>®</sup> Peel and Place Dressing Kit, Small	16.9 cm x 20.6 cm	6.1 cm x 8.6 cm	2 cm or less
EZ5MED - Solventum™ V.A.C. <sup>®</sup> Peel and Place Dressing Kit, Medium	23.7 cm x 29.2 cm	11.1 cm x 16.6 cm	4 cm or less
EZ5LRG - Solventum™ V.A.C. <sup>®</sup> Peel and Place Dressing Kit, Large	26 cm x 35.6 cm	13.6 cm x 24.2 cm	6 cm or less

## Symbols Glossary

Symbol Title	Symbol	Description and Reference
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and / or 2014/30/EU
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Container quantity		Indicates quantity of medical devices in container
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Do not resterilize		Indicates a medical device that is not to be resterilized. Source: ISO 15223, 5.2.6
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4

MR Safe		An item that poses no known hazards resulting from exposure to any MR environment. MR Safe items are composed of materials that are electrically nonconductive, nonmetallic, and nonmagnetic. Source: ASTM F2503-20 Fig. 4
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Refer to instruction manual / booklet	(*) 	To signify that the instruction manual / booklet must be read. ISO 7010-M002
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
Skin Performance		The feather designates products that deliver the securement power you need with less damage to skin.
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Sterilized using irradiation		Indicates a medical device that has been sterilized using irradiation. Source: ISO 15223, 5.2.4
Swiss Authorized Representative		Indicates the authorized representative in Switzerland. Source: Swissmedic.ch

\* This symbol is blue on packaging.

Unique device identifier	<b>UDI</b>	Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Up to 7 Days Use		Dressings may be left in place for up to 7 days.

For more information see, [elFU.Solventum.com](http://elFU.Solventum.com)

© Solventum 2025. Solventum, the S logo and V.A.C. Peel and Place are trademarks of Solventum or its affiliates.

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KCI Medical GmbH
KCI Medical Limited
KCI USA, Inc.