Acelity™

PREVENA™ INCISION MANAGEMENT SYSTEM

PREVENA™ PEEL & PLACE™ DRESSING

FOR USE WITH PREVENA™ THERAPY UNIT AND
KCI V.A.C.® THERAPY UNITS

**APPLICATION INSTRUCTIONS
FOR CLINICIANS ONLY**

Rx Only

Table of Contents

Product Description and Indication for Use	5
Important Information for Users.....	5
Optimum Use Conditions.....	5
Contraindication.....	6
Warnings.....	6
Precautions	7
Prevena™ Incision Management System Site Preparation	8
Drain Tubes and Pain Management Control Devices.....	8
Prevena™ Peel & Place™ Dressing 13 cm and 20 cm.....	9
Prevena™ Peel & Place™ Dressing Kit Components	9
Dressing Application Instructions.....	10
Using the Prevena™ Dressing with the Prevena™ 125 Therapy Unit.....	11
Connecting the Prevena™ Dressing to Prevena™ 125 Therapy Units..	11
Activating the Prevena™ Therapy Unit	11
Duration of Prevena™ Therapy	12
Dressing Removal.....	12
Prevena™ 125 Therapy Unit Indicators and Alerts.....	13
Prevena™ 125 Therapy Unit VisiCheck™ Feature	14
Correcting a Leak Condition.....	14
Check Canister and Tubing Connection.....	15
Indications That a Leak Condition Has Been Corrected	15
Prevena™ Dressing Pressure Indicator	16
Prevena™ 125 Therapy Unit Disposal	16
Battery Replacement on the Prevena™ 125 Therapy Unit	16
Instructions For Patient	17
Daily Use	17
Sleeping	17
Showering and Bathing	17
Strenuous Activity.....	17

Cleaning.....	18
Device Disposal.....	18
Using the Prevena™ Dressing with KCI V.A.C.® Therapy Units	18
Alarm Resolutions	18
Connecting the Prevena™ Dressing to V.A.C.® Therapy Units.....	18
Setting Negative Pressure on the V.A.C.® Therapy Units	18
Specifications	19
Electromagnetic Compatibility	19
Symbols Used	23
Customer Contact Information	23
Bibliography of Published Studies	24

INSTRUCTIONS FOR USE

PREVENA™ INCISION MANAGEMENT SYSTEM

KCI customer contact information is located in the back of this guide.

Product Description and Indication for Use

The Prevena™ Incision Management System is intended to manage the environment of surgical incisions and surrounding intact skin in patients at risk for developing post-operative complications, such as infection, by maintaining a closed environment via the application of a negative pressure wound therapy system to the incision. The Prevena™ Dressing skin interface layer with silver reduces microbial colonization in the fabric.

It consists of:

A Prevena™ Peel & Place™ Dressing and a source of negative pressure, which may be one of the following KCI therapy units:

- Prevena™ 125 Therapy Unit
- ActiV.A.C.® Therapy Unit
- InfoV.A.C.® Therapy Unit
- V.A.C.ATS® Therapy Unit
- V.A.C. Freedom® Therapy Unit
- V.A.C.Ultra™ Therapy Unit

Important Information for Users

As with any prescription medical device, failure to carefully read and follow all instructions and safety information prior to use may lead to improper product performance.

All of the Prevena™ Incision Management System components are disposable and are for single use only. Re-use of disposable components may result in wound contamination, infection and / or failure of the wound to heal.

CAUTION: Prevena™ Therapy should be applied and removed only by qualified physicians, nurses or caregivers.

Optimum Use Conditions

For maximum benefit the Prevena™ System should be applied immediately post surgery to clean surgically closed wounds. It is to be continuously applied for a minimum of two days up to a maximum of seven days. It can transition home with the patient; however, all Prevena™ Dressing changes should be performed under direct medical supervision.

Prevena™ Therapy will not be effective in addressing complications associated with the following:

- ischemia to the incision or incision area
- untreated or inadequately treated infection
- inadequate hemostasis of the incision
- cellulitis of the incision area

Prevena™ Therapy should not be used to treat open or dehisced surgical wounds. The V.A.C.® Therapy System should be considered for treatment of these wounds.

Prevena™ Therapy should be used with caution in the following patients:

- patients with fragile skin surrounding the incision as they may experience skin or tissue damage upon removal of the Prevena™ Dressing
- patients who are at an increased risk of bleeding from the incision associated with the use of anticoagulants and / or platelet aggregation inhibitors

Contraindication

- sensitivity to silver

Warnings

Bleeding: Before applying the Prevena™ Incision Management System to patients who are at risk of bleeding complications due to the operative procedure or concomitant therapies and / or co-morbidities, ensure that hemostasis has been achieved and all tissue planes have been approximated. If active bleeding develops suddenly or in large amounts during therapy, or if frank blood is seen in the tubing or in the canister, the patient should leave the Prevena™ Dressing in place, turn off the Prevena™ 125 Therapy Unit and seek immediate emergency medical assistance.

Infected Wounds: As with any wound treatment, clinicians and patients / caregivers should frequently monitor the patient's wound, periwound tissue and exudate for signs of infection or other complications. Some signs of infection are fever, tenderness, redness, swelling, itching, rash, increased warmth in the wound or periwound area, purulent discharge or strong odor. Infection can be serious, and can lead to complications such as pain, discomfort, fever, gangrene, toxic shock, septic shock and / or fatal injury. Some signs or complications of systemic infection are nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucus membranes, disorientation, high fever, refractory and / or orthostatic hypotension or erythroderma (a sunburn-like rash). Silver in the interface layer of the Prevena™ Dressing is not intended to treat infection, but to reduce bacterial colonization in the fabric. **If infection develops, Prevena™ Therapy should be discontinued until the infection is treated.**

Allergic Response: The Prevena™ Dressing has an acrylic adhesive coating and a skin interface layer with silver, which may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives or silver. If a patient has a known allergy or hypersensitivity to these materials, do not use Prevena™ Dressings. If any signs of allergic reaction, irritation or hypersensitivity develop, such as redness, swelling, rash, urticaria or significant pruritus, patient should consult a physician immediately. If bronchospasm or more serious signs of allergic reaction appear, the patient should turn off the therapy unit and seek immediate emergency medical assistance.

Defibrillation: Remove the Prevena™ Dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and / or patient resuscitation.

Magnetic Resonance Imaging (MRI): The Prevena™ 125 Therapy Unit is MR unsafe. Do not take the Prevena™ 125 Therapy Unit into the MR environment. The Prevena™ Dressings can typically remain on the patient with minimal risk in an MR environment. Interruption of Prevena™ Therapy during MRI may reduce the effectiveness of Prevena™ Therapy. Prevena™ Dressings pose no known hazards in an MR environment with the following conditions of use: static magnetic field of 3 Tesla or less, spatial gradient field of 720 Gauss / cm or less, and maximum whole-body-averaged specific absorption rate (SAR) of 3W / kg for 15 minutes of scanning.

Diagnostic Imaging: The Prevena™ Dressing contains metallic silver that may impair visualization with certain imaging modalities.

Hyperbaric Oxygen Therapy (HBO): Do not take the Prevena™ 125 Therapy Unit or Prevena™ Dressings into a hyperbaric oxygen chamber. They are not designed for this environment and **should be considered a fire hazard.** If Prevena™ Therapy is reinitiated after HBO treatment, do not readhere the same dressing; a new dressing must be applied.

Canister Full: If at any time while using the Prevena™ System the canister becomes full of fluid, indicated by a Maximum Capacity Alert or visual inspection, the patient should turn off the therapy unit and contact the treating physician.

Standard Operation: Do not use accessories or materials not provided with the Prevena™ Incision Management System. For a list of acceptable therapy units with which Prevena™ Dressings may be used, see the **Using the Prevena™ Dressing with KCI V.A.C.® Therapy Units** section on page 18.

Precautions

Standard Precautions: To reduce the risk of transmission of bloodborne pathogens, apply standard precautions for infection control with all patients, per institutional protocol, regardless of their diagnosis or presumed infection status.

Circumferential Dressing Application: Avoid applying the Prevena™ Dressing circumferentially. In cases where the clinician determines that the benefits of applying the Prevena™ Dressing circumferentially outweigh the risk of circulatory compromise, extreme care should be taken not to stretch or pull the dressing when securing it. Attach the dressing loosely and stabilize edges with an elastic wrap if necessary. It is crucial to systematically and recurrently palpate distal pulses and assess distal circulatory status. If circulatory compromise is suspected, discontinue therapy and remove dressing.

Electrodes or Conductive Gel: Do not allow the Prevena™ Dressing to come in contact with EKG or other electrodes or conductive gels during electronic monitoring or when taking electronic measurements.

Dressing Components: The Prevena™ System contains ionic silver (0.019%). Application of products containing silver may cause temporary tissue discoloration.

- Always use Prevena™ Dressings and canisters from sterile packages that have not been opened or damaged.
- All components of the Prevena™ Therapy are for single use only. Do not re-use any component of this system.
- To avoid trauma to the skin, do not pull or stretch the adhesive border of the dressing during application.

Prevena™ Incision Management System Site Preparation

1. Prior to surgery, shave or clip the surgical area where the dressing will be applied to improve dressing adhesion and seal integrity.
2. Gather all items needed for application:
 - sterile wound cleaning solution, e.g. water, saline or alcohol
 - sterile gauze or other material to clean application site
 - Prevena™ System (Check expiration date on box)
3. After surgery, cleanse the application site with sterile gauze and sterile wound cleaning solution using a circular motion beginning at the center of the surgical area and extending outward to ensure that the site is free of foreign material.
4. Pat the application site dry with sterile gauze. To ensure proper adhesion, the application site must be completely dry before dressing is applied.

Drain Tubes and Pain Management Control Devices


The Prevena™ System can be used with both drain tubes and pain devices, provided the dressing is not placed over tubing where it exits the skin. Surgical drains must be routed under the skin beyond the boundary of the dressing and function independently of the Prevena™ System.

NOTE: While the concomitant use of surgical drains is allowable with the Prevena™ Therapy, the system must not be used as an outlet or reservoir for the drain.


Prevena™ Peel & Place™ Dressing 13 cm and 20 cm

Prevena™ Peel & Place™ Dressing Kit Components

The Peel & Place™ Dressing Kit contains the following single use, disposable components.

Prevena™ Peel & Place™ Dressing - 13 cm with Pressure Indicator - a specially designed dressing for application to the surgical area (**13 cm kit only**)

Prevena™ Peel & Place™ Dressing - 20 cm with Pressure Indicator - a specially designed dressing for application to the surgical area (**20 cm kit only**)


Prevena™ Patch Strips™ - used to help seal leaks around dressing


V.A.C.® Connector - used to connect the Prevena™ Dressing to a KCI V.A.C.® Therapy Unit

Ruler - the removeable label may be used as needed to record the date of dressing application or removal

Prevena™ System 45 mL Canister - a sterile reservoir for collection of wound fluids


Prevena™ 125 Therapy Unit - delivers negative pressure to the surgical area. The unit is battery powered. The non-sterile Prevena™ 125 Therapy Unit Carrying Case is provided to facilitate patient mobility.

Dressing Application Instructions

(Illustrations in the steps below show Prevena™ Peel & Place™ Dressing - 20 cm)


CAUTION: If the dressing covers the umbilicus, the umbilicus must first be fully filled with an anti-microbial petroleum gauze prior to dressing application.

1. Open the sterile dressing package and remove dressing and patch strips using aseptic technique. Do not use if package has been torn or the sterile seal has been compromised.


2. Gently peel back the center strip on the back of the dressing exposing the pull tabs and adhesive.


3. Center and apply the dressing over the closed wound or incision ensuring that the adhesive will not contact or cover the surgical closure. Orient the dressing on the patient to eliminate sharp bends or kinks in the tubing.


4. Remove the remaining bottom adhesive covers by grasping the bottom tabs and gently pulling.

5. Firmly press around the dressing to ensure a good seal where the adhesive contacts the skin.

6. Remove top stabilization layers.


NOTE: The removeable label on the supplied ruler may be used as needed to record date of dressing application or removal.


Using the Prevena™ Dressing with the Prevena™ 125 Therapy Unit

Connecting the Prevena™ Dressing to Prevena™ 125 Therapy Units


1. Remove the Prevena™ System 45 mL Canister from the sterile package. Do not use if package has been torn or the sterile seal has been compromised.
2. Insert the canister into the Prevena™ 125 Therapy Unit and slide inward until canister clicks. Canister is fully inserted when the side tabs are flush with the body of the therapy unit.
3. Connect the dressing tubing to the canister tubing by twisting the connectors until they lock.
4. Begin therapy.

Activating the Prevena™ Therapy Unit

Green Light


On / Off Button


1. Press and hold the On / Off button for two seconds; an audible beep will confirm that therapy is on. A green light on the front of the unit indicates that therapy is on.

NOTE: Pressing the On / Off button begins the 192 hour (eight day expected service life) life cycle of the therapy unit. Turning the therapy unit off stops the life cycle counter. Turning the therapy unit on for purposes other than delivering therapy reduces the life cycle of the therapy unit. It is not recommended to press the On / Off button until therapy is ready to begin.


NOTE: To turn therapy unit off, press and hold the On / Off button for five seconds.

2. With therapy on, assess dressing to ensure integrity of seal.
 - The dressing should have a wrinkled appearance and the foam bolster should be compressed.
 - The pressure indicator on the dressing should be in the collapsed position.
 - If there is any evidence of a leak, refer to the **Prevena™ 125 Therapy Unit Indicators and Alerts** section on page 13.

3. Place the therapy unit into the Prevena™ 125 Therapy Unit Carrying Case. Make sure the display is visible through the opening in the carrying case when the front flap is lifted.

4. The Prevena™ 125 Therapy Unit Carrying Case has an integrated belt loop and a separate adjustable strap to allow for versatile positioning.

CAUTION: Do not wrap carrying case strap or dressing tubing around neck.



Duration of Prevena™ Therapy

- Therapy should be continuous for a minimum of two days up to a maximum of seven days.
- The Prevena™ 125 Therapy Unit will automatically time-out after 192 hours (eight days) of cumulative run time.
- Patients should be instructed to contact their treating physician and not to turn therapy off unless:
 - advised by the treating physician
 - bleeding develops suddenly or in large amounts during therapy
 - there are signs of allergic reaction or infection
 - the canister is full of fluid
 - batteries need to be changed
 - system alerts must be addressed
- Patient should be instructed to contact the treating physician if therapy unit turns off and cannot be restarted before therapy is scheduled to end, or if canister becomes full of fluid.
- At end of therapy, patient should return to treating physician for dressing removal.

Dressing Removal

NOTE: If dressing is lifted to observe wound, do not re-adhere the same dressing; a new dressing must be applied.

WARNING: Dressings should always be removed in-line with the sutures and NEVER across the sutures.

1. Turn Prevena™ 125 Therapy Unit off by pressing and holding the On / Off button for five seconds.
2. Gently stretch the drape / dressing horizontally to release the adhesive from the skin. Do not peel vertically. Remove the drape / dressing in-line with the sutures, NEVER across the sutures.
3. Clean any residual adhesive with an alcohol swab.

If a new dressing is to be applied:

1. Ensure the incision area is clean by using an alcohol swab or an antiseptic wipe.
2. Allow skin to dry completely before applying a new drape.
3. Follow the **Dressing Application Instructions** on page 10.

Prevena™ 125 Therapy Unit Indicators and Alerts


Caution Indicator

Audible Alarm (Muted)

Battery Level

Visual Alerts - Solid LEDs cannot be turned off by the user. Visual alerts will only stop when the alert condition has been corrected.

Audible Alerts - Repeated beeps (which in some cases will increase in volume) can be temporarily muted (paused) by pressing the On / Off button once. The audible alert will re-occur after 60 minutes unless the alert condition has been corrected.

The Prevena™ 125 Therapy Unit will provide audible and visual alerts as shown below.

If alert conditions cannot be corrected, patient should contact treating physician. For further product support see **Customer Contact Information** section on page 23.

Leak Alert

One beep, one solid yellow light.

See **Correcting a Leak Condition** section on page 14.

When a leak condition is corrected the therapy unit will cancel alert. There may be a delay between when the leak is corrected and when the alert is discontinued.

Canister Full Alert

Two beeps, one solid yellow light.

Visually inspect canister. If full or near full, turn therapy unit off and call the treating physician immediately.

Low Battery Alerts

Low Battery Alert - One slow beep, one solid yellow light. Change batteries within six hours.

Critical Battery Alert - One beep, repeating rapidly, increasing in volume, one solid yellow light. Change batteries immediately.

Battery replacement will cancel the alert. See the **Battery Replacement on the Prevena™ 125 Therapy Unit** section on page 16.

System Error Alert

One beep, repeating rapidly, increasing in volume, two solid yellow lights. Turn the unit off and then on again.

If alert continues, notify the treating physician.

Device Lifecycle Expired

Three solid yellow lights. If the lifecycle expires while therapy unit is on, a beep will sound for 15 seconds and then automatically shut off.

If the therapy unit is off, has exceeded the run time, and an attempt is made to turn the therapy unit on, the therapy unit will alert for three seconds and automatically shut off. Call the treating physician.

Prevena™ 125 Therapy Unit VisiCheck™ Feature


To ensure proper application of the Prevena™ Dressing, the Prevena™ 125 Therapy Unit provides the VisiCheck™ Feature.

By double pressing the On / Off button, the unit will display the leak rate of the system for three seconds.

To prevent nuisance leak alarms, the leak rate status should be "Best" (one light illuminated) or "Good" (two lights illuminated).


If the VisiCheck™ Feature indicates a "Marginal" (three lights illuminated) leak rate condition, refer to the following **Correcting a Leak Condition** section below for methods of reducing the leak rate of the system. The system leak rate is calculated every seven seconds. If a corrective action is taken to reduce a leak rate, use the VisiCheck™ Feature afterward to verify the leak rate condition was corrected.

When the therapy unit detects a leak:


- A visual and audible leak alert will activate. See the **Prevena™ 125 Therapy Unit Indicators and Alerts** section on page 13.
- Therapy unit will turn on more frequently.

Correcting a Leak Condition

(Illustrations in the steps below show Prevena™ Peel & Place™ Dressing - 20 cm)


1. With therapy unit on, slowly press firmly around dressing edge to ensure good contact between adhesive and skin.


2. If a leak is identified, use Prevena™ Patch Strips™ (located in dressing package) to help seal leaks around dressing. If large wrinkles are present, place patch strips so they run in line along the length of the wrinkle and not across the wrinkle.

Check Canister and Tubing Connection

1. Ensure canister is securely locked in the therapy unit. When canister is installed, a distinct click will be heard indicating it has been properly installed. Side tabs on canister should be flush with unit.


2. Check tubing connectors to ensure they are fully engaged and locked.


Indications That a Leak Condition Has Been Corrected

- Therapy unit will become quiet, running only intermittently.
- Audible leak alert will stop; visual alert will turn off. There may be a brief delay between when the leak is corrected and the alert is discontinued.


Prevena™ Dressing Pressure Indicator

Each dressing is equipped with a pressure indicator that shows when system pressure is at an acceptable level.

Indicator up - system pressure not acceptable.

Indicator down - system pressure acceptable.

Prevena™ 125 Therapy Unit Disposal


At the end of therapy, the patient should return the Prevena™ 125 Therapy Unit to the physician for disposal. Dispose of all waste according to local requirements. Improper disposal may run the risk of regulatory non-compliance.

Battery Replacement on the Prevena™ 125 Therapy Unit


Battery replacement should be done as quickly as possible after a low battery alert to prevent therapy down time.

1. Turn therapy unit off (press and hold On / Off button for five seconds).


2. Remove therapy unit from carrying case. Turn the therapy unit over to expose the back side of the unit.

3. Locate the battery door and push to slide open. Install three AA batteries (lithium batteries are recommended for optimal performance) into the battery compartment.

NOTE: Always replace with new batteries. Do not mix new batteries with used batteries.

NOTE: The inside of the battery compartment is stenciled with "+" positive and "-" negative to aid in proper battery installation.

4. Close battery door.
5. Return therapy unit to carrying case.
6. Turn therapy unit on to resume therapy (press and hold On / Off button for two seconds).

Instructions For Patient

Review the following information with the patient prior to discharge. This information is summarized in the Prevena™ Therapy System Patient Guide which must be provided to the patient at discharge.

Daily Use

The Prevena™ Therapy is portable and small enough that it may be worn beneath clothing during normal patient activities as approved by the treating physician.

CAUTION: Advise patient to NOT SUBMERGE therapy unit or dressing in liquid and to ensure therapy unit is not pulled into a tub or sink where it may become submerged.

CAUTION: The Prevena™ 125 Therapy Unit is a medical device not a toy. Keep away from children, pets and pests as they can damage the dressing and therapy unit and affect performance. Keep therapy unit free of dust and lint.

Sleeping

Instruct patient to:

- place the therapy unit in a position where tubing will not become kinked or pinched.
- ensure therapy unit will not be pulled off a table or fall to the floor during sleep.

Showering and Bathing

Advise patient of the following recommendations:

- Light showering is permissible, bathing is not.
- While showering, the device and dressing should be protected from prolonged direct spray and / or being submerged.
 - Hanging the therapy unit from a soap / shampoo holder or shower head may be acceptable if protected from prolonged direct spray.
 - Dressing may be exposed to common shower soaps and rinsed with indirect shower stream. Do not submerge dressing. Do not remove dressing.
- When towel drying, avoid disturbing or damaging dressing.

Strenuous Activity

Advise patient as to when and at what level physical activities may be resumed. It is recommended that patients avoid strenuous activity while using the Prevena™ Therapy System.

Cleaning

Advise patient that the Prevena™ 125 Therapy Unit and Prevena™ 125 Therapy Unit Carrying Case can be wiped with a damp cloth using a mild household soap solution that does not contain bleach.

Device Disposal

At the end of therapy, the patient should return therapy unit to physician for disposal. Dispose of all waste according to local requirements. Improper disposal may run the risk of regulatory non-compliance.

Using the Prevena™ Dressing with KCI V.A.C.® Therapy Units

In certain cases as directed by the treating physician, Prevena™ Dressings can be used with a KCI V.A.C.® Therapy Unit.

The ActiV.A.C.®, InfoV.A.C.®, V.A.C.ATS®, V.A.C.Ultra™ and V.A.C. Freedom® Therapy Units may be used on closed surgical incisions. Refer to the therapy unit's user manual for complete instructions for use.

WARNING: Review all negative pressure wound therapy system safety information before initiating therapy.

WARNING: DO NOT use with V.A.C. VeraFlo™ Therapy (Instillation) provided by the V.A.C.Ultra™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

NOTE: When Prevena™ Dressings are used with V.A.C.® Therapy, the SensaT.R.A.C.™ Dressing wound pressure sensing function on the therapy unit will be disabled. Instead, the Prevena™ Dressing Pressure Indicator (page 16) should be used to verify that incision site pressure is acceptable.

Alarm Resolutions


KCI V.A.C.® Therapy Unit alarms should be addressed in a timely manner. Refer to the appropriate therapy unit user manual for complete information on alarm resolutions.

Connecting the Prevena™ Dressing to V.A.C.® Therapy Units

V.A.C.® Connector

A KCI V.A.C.® Connector will be needed to connect the Prevena™ Dressing to a KCI V.A.C.® Therapy Unit. This connector, which is available in the Prevena™ Dressing package, **must** be used for negative pressure wound therapy to work effectively and accurately. No other method for connecting the dressing to the therapy unit should be attempted.

1. Connect the Prevena™ Dressing tubing to the V.A.C.® Therapy Unit canister tubing using the V.A.C.® Connector:

- Push the connectors together
- Twist the connectors to lock

2. Ensure clamp on canister tubing is open.

Setting Negative Pressure on the V.A.C.® Therapy Units

Activate V.A.C.® Therapy at -125mmHg continuous. Do not choose any other negative pressure setting, intermittent or DPC modes of negative pressure.

Specifications

Storage Conditions

Temperature Range - 4°F (-20°C) to 140°F (60°C)
 Relative Humidity Range 15% - 95%, non-condensing

Operating Conditions

Temperature Range 41°F (5°C) to 104°F (40°C)
 Atmospheric Pressure Range 1060 hpa (-1253 ft /-381.9 m)
 (for Optimum Performance) to 700 hpa (9878 ft / 3010 m)

Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.

The Prevena™ Therapy System is classified as a Type BF applied part under IEC 60601-1:2005 (3rd Edition).

IP24 - Protection against solid objects greater than 12.5 mm and against liquid water sprays for short periods of time.

All alerts are classified as low priority according to IEC 60601-1-8:2006.

Electromagnetic Compatibility

Electromagnetic Interference - Although this equipment conforms with the intent of the directive 89 / 336 / EEC in relation to electromagnetic compatibility (EMC), all electrical equipment may produce interference. If interference is suspected, move equipment away from sensitive devices or contact the manufacturer.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the following tables.

The following tables document compliance levels and guidance from the IEC 60601-1-2:2007 Standard, for the electromagnetic environment in which the Prevena™ 125 Therapy Unit should be used. The Prevena™ 125 Therapy Unit also meets the standard for electromagnetic compatibility related to use in the home care environment (IEC 60601-1-11:2010).

Guidance and manufacturer's declaration - electromagnetic emissions		
The Prevena™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or user of the Prevena™ 125 Therapy Unit should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic environment
RF emissions CISPR 11	Group 1	The Prevena™ 125 Therapy Unit uses RF energy only for internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. Battery operated device
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration - electromagnetic immunity

The Prevena™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or user of the Prevena™ 125 Therapy Unit should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	In accordance with IEC 60601-1-2: 2007, floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	Not Applicable	Battery operated device
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not Applicable	Battery operated device
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_r (>95 % dip in U_r) for 0.5 cycle 40 % U_r (60 % dip in U_r) for 5 cycles 70 % U_r (30 % dip in U_r) for 25 cycles <5 % U_r (>95 % dip in U_r) for 5 cycles	Not Applicable	Battery operated device
Power frequency (50Hz / 60Hz) magnetic field IEC 61000-4-8	3 A / m	3 A / m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_r is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the Prevena™ 125 Therapy Unit

The Prevena™ 125 Therapy Unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Prevena™ 125 Therapy Unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Prevena™ 125 Therapy Unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz not applicable	80 MHz to 800 MHz $d = \left[\frac{3.5}{E1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E1} \right] \sqrt{P}$
0.01	not applicable	0.12	0.23
0.1	not applicable	0.37	0.74
1	not applicable	1.17	2.33
10	not applicable	3.69	7.38
100	not applicable	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separate distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1, At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from surfaces, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity

The Prevena™ 125 Therapy Unit is intended for use in an electromagnetic environment specified below. The customer or user of the Prevena™ 125 Therapy Unit should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not Applicable	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Prevena™ 125 Therapy Unit, including cables, than the recommended separation distance calculated from the equation application to the frequency of the transmitter.</p> <p style="text-align: center;">Recommended separation distance</p> <p style="text-align: center;">Battery Operated Device</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹, should be less than the compliance level in each frequency range².</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3 V / m 80 MHz to 2.5 GHz	3 V / m 80 MHz to 2.5 GHz	

NOTE 1, At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2, These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

¹ Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Prevena™ 125 Therapy Unit is used exceeds the applicable RF compliance level above, the Prevena™ 125 Therapy Unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Prevena™ 125 Therapy Unit.

² Over the frequency range 150kHz, field strengths should be less than $[V] \text{ V / m}$.

Symbols Used

Do not use if package is damaged or open

IP24

Ingress Protection

Type BF applied part

STERILE R

Sterile using radiation

Date of Manufacture

Manufacturer

Fragile

Keep Dry

Each

Content Information

Temperature Limits

ETL Listed
Conforms to AAMI STD ES60601-1
Certified to CSA STD C22.2 No. 60601-1

Conforms with the Medical Device Directive (93/42/EEC) and has been subject to the conformity procedures laid down in the council directive

Rx only

CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.

Consult Instructions for Use

Refer to Clinician Guide

Single Use Only

EC REP

Authorized Representative in the European Community

Non Sterile

Use By

REF

Catalog Number

LOT

Lot Number

Do Not Resterilize

This product is designated for separate collection at an appropriate collection point. Do not dispose of as household waste.

Customer Contact Information

For questions regarding this product, supplies, maintenance or additional information about KCI products and services, please contact KCI or a KCI authorized representative, or:

Outside the US visit www.kci-medical.com

In the US call 1-800-275-4524 or visit www.kci1.com

KCI USA, Inc. 12930 IH10 West, San Antonio, TX 78249

Bibliography of Published Studies

- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.

- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

KCI USA, Inc.
12930 IH 10 West
San Antonio, TX
78249 USA
www.kci1.com

KCI Medical Products (UK), Ltd.
11 Nimrod Way
Wimborne, Dorset, BH21 7SH
United Kingdom
www.kci-medical.com

All trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and / or licensors.
©2015 KCI Licensing, Inc. All rights reserved. 416356 Rev A 6 / 2015

Acelity™

DE

PREVENA™ INCISION MANAGEMENT SYSTEM

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 16. Juni 2009(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. Feb. 2011;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 1. Dezember 2013;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30. Januar 2013
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 19. Feb. 2013 [Elektronische Veröffentlichung vor dem Druck]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 1. Dezember 2013;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 6. Dezember 2011;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 1. Okt. 2013;71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 21. Aug. 2012(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 1. September 2013;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 24. Nov. 2012 (11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. Nov.-Dez. 2006;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. Mai 2013;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. Jan./Feb. 2011;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 20. September 2013;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 1. Dezember 2013;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. Sep. 2011;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. Juni 2012;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. März 2013;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. Apr. 2012;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. Juni 2013;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. Juni 2010;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. Juni 2006;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. Jan. 2012;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. März 2013;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. Dez. 2012;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. März 2012;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 juni 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 februari;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 december 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 januari 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 februari 19. [Epub voorafgaand aan gedrukte publicatie]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 december 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 december 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 oktober 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 augustus 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 september 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 november 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 november-december;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 mei;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 januari/februari;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 september 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 december 1;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 september;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 2013 september 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 juni;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 maart;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 april;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 juni;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 juni;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 juni;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 januari;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 maart;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 december;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 maart;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb 8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.

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- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. Junio de 2009 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 1 de diciembre de 2013;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30 de enero de 2013.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Epub en impresión].
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 1 de diciembre de 2013;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 6 de diciembre de 2011;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 1 de septiembre de 2013;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dic;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 20 de septiembre de 2013;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 1 de diciembre de 2013;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 20 de septiembre de 2013.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. Juni 2009 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. December 2013 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30. Januar 2013.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 19. Feb. 2013. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. December 2013 1;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. December 2011 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery Okt 2013 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. Aug 2012 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. September 2013 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. Nov 2012 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. Nov-Dec 2006;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. Maj 2013;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. Jan/Feb 2011;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty September 2013 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. December 2013 1;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. Sep 2011;19(5):588-96.

- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations September 2013 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. Jun 2012;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. Mar 2013;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. Apr 2012;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. Jun 2013;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. Jun 2010;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. Jun 2006;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. Jan 2012;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. Mar 2013;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. Dec 2012;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. Mar 2012;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 Mar;16(2):140-6
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 16 de junho de 2009(2): 140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. Fev 2011;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 1 de dezembro de 2013; 12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30 de janeiro de 2013.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 19 de fevereiro de 2013. [publicado eletronicamente antes de impresso]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 1 de dezembro de 2013;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 6 de dezembro de 2011;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery Out 2013 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 21 de agosto de 2012(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 1 de setembro de 2013;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 24 de novembro de 2012(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. Nov-Dez 2006;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. Maio de 2013;145(5):1387-92.
- Haghshenas-kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. Jan/Fev 2011;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 20 de setembro de 2013;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 1 de dezembro de 2013;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. Set de 2011;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 20 de setembro de 2013.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. Jun 2012;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. Mar 2013;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. Abr 2012;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. Jun 2013;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. Jun 2010;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. Jun 2006;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. Jan 2012;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. Mar 2013;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. Dez 2012;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. Mar 2012;19(1):67-75.

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- Atkins BZ, v.d. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 16 Haziran 2009 (2):140-6.
- Atkins BZ, v.d. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Şubat;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 1 Aralık 2013;12(4):172-4.
- Blackham AU, v.d. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30 Ocak 2013.
- Bollero D, v.d. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 19 Şubat 2013. [Epub ahead of print]
- Bonds AM, v.d. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 1 Aralık 2013;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 6 Aralık 2011;6(1):160.
- Conde-Green A, v.d. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery Ekim 2013 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 21 Ağustos 2012 (8):386-8.
- Faroqi L, v.d. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 1 Eylül 2013;5(2):44-7.
- Glaser DA, v.d. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 24 Kasım 2012 (11):308-11
- Gomoll AH, v.d. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Kasım-Aralık;20(10):705-9.
- Grauhan O, v.d. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. Mayıs 2013;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. Ocak/Şubat 2011;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 20 Eylül 2013;13:e49.
- Karlakki S, v.d. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 1 Aralık 2013;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. Eylül 2011;19(5):588-96.

- Mark KS, v.d. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 20 Eylül 2013.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. Haziran 2012;255(6):1043-7.
- Matatov T, v.d. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. Mart 2013;57(3):791-5.
- Pachowsky M, v.d. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. Nisan 2012;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. Haziran 2013;14(3):270-4.
- Reddix RN, v.d. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. Haziran 2010;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. Haziran 2006;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. Ocak 2012;26(1):37-42.
- Tauber R, v.d. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. Mart 2013;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. Aralık 2012;204(6):1021-4.
- Wilkes RP, v.d. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. Mart 2012;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Julkaistu ensin sähköisessä muodossa]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faroqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghshenasskashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

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Alle andre varemerker som nevnes her, tilhører KCI Licensing, Inc., deres datterselskaper og/eller lisensgivere.
©2015 KCI Licensing, Inc. Med enerett. 416356 Rev A 6 / 2015