
# KerraCel™

## Gelling fiber dressing

# Instructions for use

#### Information for use

#### Product Description

KerraCel™ is a sterile, gelling fiber dressing, made from carboxymethyl cellulose (CMC). The non-woven gelling fibers absorb exudate from the wound, turning them into a soft gel which maintains a moist wound healing environment.

#### Indications

Suitable for use on moderate and heavily exuding wounds (chronic or acute) such as:

- Leg ulcers
- Pressure ulcers
- Diabetic foot ulcers
- Surgical wounds
- Donor sites
- Cavity wounds
- Superficial and partial thickness burns
- Trauma wounds
- Oncology wounds
- Control minor bleeding in superficial wounds

KerraCel is suitable for use under compression bandaging.

#### Contra-indications

Known allergy or hypersensitivity to the ingredient of KerraCel, carboxymethyl cellulose. KerraCel should not be used:

- For surgical implantation
- To control heavy bleeding – the dressing is not a haemostat

#### Precautions/Warnings

- For single use only
- Do not re-use. Re-using will cause cross contamination
- In case of clinical infection, consult a healthcare professional
- Dispose in normal waste is used in a household
- Dispose in clinical waste if used in a clinical establishment

#### Using KerraCel

#### Exuding Wounds

- Cleanse the wound in accordance with local procedures.
- Select an appropriate dressing size, so the dressing overlaps the wound margin by approximately 1cm/0.4”
- Open the pack and lay the sterile dressing onto the wound. Either side of the sterile dressing may be used.
- Apply the dressing to the wound and cover with an appropriate secondary dressing.
- When using the ribbon in cavity wounds, leave at least 2.5cm/1” outside the wound for ease of removal. Only pack cavity wounds by 80%, as the dressing will expand to fill the space on contact with exudate.

#### Dry Wounds

- Cleanse the wound in accordance with local procedures.
- Select an appropriate dressing size, so the dressing overlaps the wound margin by approximately 1cm/0.4.”
- Open the pack and lay the sterile dressing onto the wound. Either side of the sterile dressing may be used. Wet the dressing with saline over the wound area only.
- Cover with a suitable secondary dressing that retains the added moisture (occlusive dressing).

#### Partial Thickness Burns

- Cleanse the wound in accordance with local procedures.
- Select an appropriate dressing size, so the dressing overlaps the wound margin by approximately 5cm/2” as the dressing will shrink in size as it absorbs the excess exudate.
- Open the pack and lay the sterile dressing onto the wound. Either side of the sterile dressing may be used.
- Apply the dressing to the wound and cover with an appropriate secondary dressing.

#### Frequency of Change

Check the dressing at regular intervals and change according to the clinical condition of the wound. If the condition of the wound deteriorates, discontinue use and consult your healthcare professional. KerraCel can be left in place for several days. However, this depends upon the condition of the wound and the skin in the peri-wound area. At dressing change, if KerraCel does not come away from the wound easily, apply saline to the dressing and gently remove.

#### Notes

KerraCel can be ordered from authorised medical distributors and from pharmacies.

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