

SPACE REPRESENTS
OVERPRINT/OVELABELLING
LABEL DIMENSIONS
ILLUSTRATED ARE
25mm (H) X 55mm (W)

PRODUCT: TOPPER™ 12

CODE: TS1102.

PRODUCT DESCRIPTION: Swabs.

LANGUAGES: English, German French, Italian, Spanish, Swedish, Dutch, Portuguese, Finish, Danish & Greek.

SIZE: 140mm (W) x 120mm (D) x 260mm (H).

COLOURS: Black.

DATE: 11.10.2017

AVC #: 21102-0

BARCODE(S): Primary & Secondary.

SYMBOLGY: GS1 128 & Code 39 Mod 11.

VERSION NUMBER: 07.

MASTER

11.10.2017

Black

21102-0

EN SYMBOLS USED ON LABELLING
FR SYMBOLES UTILISÉS SUR LES ÉTIQUETTES
ES SÍMBOLOS EMPLEADOS EN EL ENVASE
PT SÍMBOLOS UTILIZADOS NOS RÓTULOS



Do not use if package is damaged.
Ne pas utiliser si l'emballage est endommagé.
No usar si el envase está dañado.
Não utilizar se a embalagem estiver danificada.



Do not resterilize.
Ne pas restériliser.
No volver a esterilizar.
Não reesterilizar.

STERILEEO

Sterilized using ethylene oxide.
Stérilisé à l'oxyde d'éthylène.
Esterilizado con óxido de etileno.
Esterilizado com óxido de etileno.

LOT

Batch number.
Numéro du lot.
Número de lote.
Número de lote.



Do not re-use (single use).
Ne pas réutiliser (à usage unique).
No reutilizar (para un solo uso).
Não reutilizar (utilização única).

REF

Catalogue number.
Numéro de catalogue.
Número de catálogo.
Número de catálogo.

CE
0123

CE-mark and identification number of notified body.
Marquage CE et numéro d'identification de l'organisme accrédité.
Marcado CE y Número de identificación del organismo notificado.
Marcação CE e número de identificação do organismo notificado.



Use by: year, month, day.
Utiliser avant le : année, mois, jour.
Fecha de caducidad: año, mes, día.
Utilizar até: ano, mês, dia.

EC REP

Authorized representative in the European Community.
Représentant agréé dans la Communauté européenne.
Representante autorizado en la Comunidad Europea.
Representante autorizado na Comunidade Europeia.



Manufacturer.
Fabricant.
Fabricante.
Fabricante.

TOPPER™ 12

Swabs
Kompresen
Compresses
Medicazioni
Compresas
Kompreser
Compressas
Taitokset
Γάζες

Acelity™

STERILEEO

10cm x 10cm
70x2

REF TS1102

TOPPER™ 12

Swabs
Kompresen
Compresses
Medicazioni
Compresas
Kompreser
Compressas
Taitokset
Γάζες

Acelity™

STERILEEO

10cm x 10cm
70x2

REF TS1102

10cm x 10cm
70x2

STERILEEO

REF TS1102

Acelity™

Swabs
Kompresen
Compresses
Taitokset
Compresas
Medicazioni
Compresas
Γάζες


TOPPER™ 12



PZN - 08821305



(01)25051978001819



Distributed by/Distribué par/Distribuido por/Distribuido por:
©2017 Systagenix Wound Management Limited,
Gargrave, North Yorkshire, BD23 3RX, UK

EC REP Shanghai International Holding Corp. GmbH (Europe),
Eiffestrasse 80, 20537 Hamburg GERMANY

Allmed Medical Products Co., Ltd.
No.180 Gong Yuan Road, Majiadian Town,
443200 Zhijiang City, Hubei Province,
PEOPLE'S REPUBLIC OF CHINA

Systagenix
An Acelity Company



CE
0123

